## REFERENCE INTERVALS<sup>(5)</sup>

Newborn: 7-40 U/L (37°C)  
Adult: 10-35 U/L (37°C)

ALT activity in men are generally higher than levels in women.

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi<sup>®</sup> analyzer unless otherwise stated.

## RESULTS

ALT activity is reported as U/L.

## REPORTABLE RANGE (CLSI EP6)<sup>(7)</sup>

The linearity of the procedure described is 600 U/L. The lower limit of detection of the procedure described is 10 U/L. This data results in a reportable range of 10 to 600 U/L.

## PRECISION STUDIES (CLSI EP5)<sup>(7)</sup>

Data was collected on two levels of a control sera using a single lot of reagent in 40 runs conducted over 20 days.

Concentration U/L	Total SD U/L	Total CV %	Within Run SD U/L	Within Run CV %
46	1.7	3.6	1.1	2.4
144	3.5	2.5	1.1	0.7

## ACCURACY (CLSI EP9)<sup>(7)</sup>

The performance of this method (y) was compared with the performance of a similar ALT method (x) on a Roche/Hitachi<sup>®</sup> analyzer. Forty-nine patient serum samples ranging from 12 to 304 U/L were tested and gave a correlation coefficient of 0.9941. Linear regression analysis gave the following equation:

$$\text{This method} = 1.06 (\text{reference method}) - 2.8 \text{ U/L.}$$

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**SEKISUI**  
DIAGNOSTICS

**The Americas**  
Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265  
Fax: 902-628-6504  
Email: questions@sekisuidiagnostics.com  
peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

**International**  
Sekisui Diagnostics (UK) Limited  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, UK

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Temperature limitation  
Límites de temperatura

## REFERENCES/ REFERENCIAS

1. Henley, K.S., Pollard, H.M. *A New Method for the Determination of Glutamic Oxalacetic and Glutamic Pyruvic Transaminase in Plasma*, J. Lab. and Clin. Med., 46, 785 (1955).
2. Wroblewski, F., Ladue, J.S. *Serum Glutamic Pyruvic Transaminase in Cardiac and Hepatic Disease*, Proc. Soc. Exp. Biol. Med. 91, 569 (1956).
3. Bergmeyer, H.U., and Horder, M. *IFCC Methods for the Measurement of Catalytic Concentration of Enzymes, Part 3 IFCC Method For Alanine Aminotransferase*, J. Clin. Chem. Clin. Biochem. 18, 521 (1980).
4. Lines, J.G., Dennis, P.M., deCediel, N., Khayat, M.H., Krawczynski, M.J., (Editors), *AST and ALT Methods Approved*, IFCC News, No. 40 (1985/2), p. 10.
5. Burtis, C.A., Ashwood, E.R. (Ed.), *Tietz Textbook of Clinical Chemistry*, W.B. Saunders Co., Toronto, p. 790-791, 797 (1994).
6. Kaplan, Lawrence A. and Pesce, Amadeo J., (Ed), *Clinical Chemistry: Theory, Analysis, Correlation, 3<sup>rd</sup> ed.* Mosby-Year Book Inc., St. Louis, Missouri (1996), p. 518.
7. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
8. Young, D.S. *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, Third Edition (1990).

Authorized Representative/ Representante autorizado:  
Sekisui Diagnostics (UK) Limited  
50 Gibson Drive  
Kings Hill, West Malling  
Kent, ME19 4AF  
United Kingdom/ Reino Unido  
Tel (+44)(0)1732-220022  
Fax (+44)(0)1732-220024

IN31810-10  
November 3, 2011