EN

## ASPARTATE AMINOTRANSFERASE-SL ASSAY

**CATALOGUE NUMBER:** 319-10    **SIZE:** R1: 1 x 100 mL, R2: 1 x 25 mL  
319-30                            R1: 3 x 100 mL, R2: 1 x 75 mL  
319-50A                         R1: 1 x 1000 mL  
319-50B                         R2: 1 x 300 mL

### INTENDED USE


For the IN VITRO quantitative measurement of aspartate aminotransferase (AST) in serum.

### TEST SUMMARY

Karmen<sup>(1)</sup> introduced a method for the measurement of aspartate aminotransferase, formerly called glutamate oxalacetic transaminase (GOT). The procedure used in this assay incorporates several significant modifications including optimization of substrate concentrations and replacement of phosphate with tris as buffer. These modifications meet the recommendations outlined by the International Federation of Clinical Chemistry (IFCC) for AST determinations.<sup>(2), (3)</sup>

Elevated concentrations of aspartate aminotransferase result from a myocardial infarction, or hepatic disease. In cases of hepatic necrosis, elevations of AST concentrations occur prior to the onset of clinical symptoms such as jaundice. To a lesser extent elevated concentrations are indicative of diseases of some internal organs.<sup>(4)</sup>

### TEST PRINCIPLE

The AST catalyzes the reaction of alpha-ketoglutarate and L-aspartate to L-glutamate and oxalacetate. In reaction II, malate dehydrogenase (MDH) catalyzes the oxidation of NADH to NAD. The rate of decrease in absorbance of the reaction mixture at 340 nm, due to the oxidation of NADH is directly proportional to the AST activity.

### REAGENTS

AST-SL Enzyme Reagent (R1): A solution containing a buffer (pH 7.8 at 25°C), 240 mmol/L L-aspartate, >600 U/L MDH (microbial), > 600 U/L lactate dehydrogenase (microbial), and a preservative.

AST-SL Substrate Reagent (R2): A solution containing 12 mmol/L 2-oxoglutarate, > 0.18 mmol/L NADH, and a preservative.

### WARNING AND PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.  
See Material Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use. A single working reagent may be prepared by mixing 4 parts R1 with 1 part R2.

The reagents included are stable until the expiry date stated on the vial label at 2-8°C. The working reagent is stable for 4 weeks at 2-8°C.

Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum.

### SAMPLE STORAGE

AST is stable for 3 days at room temperature (approximately 25°C), 1 week at 0-4°C, and 1 month at -25°C.<sup>(5)</sup>

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(8)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

This procedure measures total AST. Red blood cells contain high concentrations of AST; therefore, hemolysis can elevate results.

Interferences from icterus and lipemia were evaluated for this method on a Roche/Hitachi<sup>®</sup> analyzer using a significance criterion of >10% variance from control.

Concentration of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant	
63 U/L	Bilirubin	40 mg/dL	684 µmol/L
112 U/L	Intralipid	400 mg/dL	1200 mg/dL (13.6 mmol/L) Simulated Triglycerides

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(6)</sup>

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Sekisui Diagnostics' Aspartate Aminotransferase (AST-SL) reagents.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material. (If applicable.)
3. Quality Control materials.

#### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:14:3.5 and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

The frequency of calibration, if necessary, using an automated system is dependent on the system and the parameters used.

#### QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

#### CALCULATIONS

The analyzer automatically calculates the aspartate aminotransferase concentration of each sample.

#### TEST LIMITATIONS

A sample with a aspartate aminotransferase concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

#### REFERENCE INTERVALS<sup>(7)</sup>

5-34 U/L (37°C)

Newborn and infant values are approximately two times higher than adult values.

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® analyzer unless otherwise stated.

## RESULTS

Aspartate aminotransferase activity is reported as U/L.

## REPORTABLE RANGE (CLSI EP6)<sup>(8)</sup>

The linearity of the procedure described is 600 U/L. The lower limit of detection of the procedure described is 10 U/L. This data results in a reportable range of 10-600 U/L.

## PRECISION STUDIES (CLSI EP5)<sup>(8)</sup>

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration U/L	Total SD U/L	Total CV %	Within Run SD U/L	Within Run CV %
54	1.7	3.2	1.2	2.1
225	5.6	2.5	1.1	0.7

## ACCURACY (CLSI EP9)<sup>(8)</sup>

The performance of this method (y) was compared with the performance of a similar method (x) on a Roche/Hitachi® analyzer. Fifty-five patient serum samples ranging from 13-301 U/L gave a correlation coefficient of 0.9979. Linear regression analysis gave the following equation:

$$\text{This method} = 1.02 (\text{reference method}) - 2.6 \text{ U/L.}$$

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**The Americas**  
Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265  
Fax: 902-628-6504  
Email: questions@sekisuidiagnostics.com  
peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

**International**  
Sekisui Diagnostics (UK) Limited  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, UK

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Temperature limitation  
Límites de temperatura

## REFERENCES/ REFERENCIAS

1. Karmen, A., Wroblewski, F., Ladue, J. Transaminase Activity in Human Blood, *J. Clin. Invest.* 34, 126 (1955).
2. Bergmeyer, H.U., Bowers, G.N., Horder, M., Moss, D.W. Provisional Recommendations on IFCC Methods for the Measurement of Catalytic Concentrations of Enzymes, Part 2 IFCC Method For Aspartate Aminotransferase, *Clin. Chem.* 23, 887 (1977).
3. Lines, J.G., Dennis, P.M., deCediel, N., Khayat, M.H., Krawczynski, M.J., (Editors), *AST and ALT Methods Approved*, IFCC News, No. 40 (1985/2), p. 10.
4. Burtis, C.A., Ashwood, E.R. (Ed.), *Tietz Textbook of Clinical Chemistry*, W.B. Saunders Co., Toronto, pp 790-791 (1994).
5. Kaplan, Lawrence A. and Pesce, Amadeo J., (Ed), *Clinical Chemistry: Theory, Analysis, Correlation*, 3<sup>rd</sup> ed. Mosby-Year Book Inc., St. Louis, Missouri (1996), p. 1072.
6. Young, D.S. *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, Third Edition (1990).
7. Kaplan, Lawrence A. and Pesce, Amadeo J., (Ed), *Clinical Chemistry: Theory, Analysis, Correlation*, 3<sup>rd</sup> ed. Mosby-Year Book Inc., St. Louis, Missouri (1996), p. 523.
8. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.

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Sekisui Diagnostics (UK) Limited  
50 Gibson Drive  
Kings Hill, West Malling  
Kent, ME19 4AF  
United Kingdom/ Reino Unido  
Tel (+44)(0)1732-220022  
Fax (+44)(0)1732-220024

IN31910-10  
July 4, 2011