EN

ALANINE AMINOTRANSFERASE L3K[®] ASSAY

CATALOGUE NUMBER: 343-25

SIZE: 2 x 125 mL

INTENDED USE



For the IN VITRO quantitative measurement of alanine aminotransferase activity in serum.

TEST SUMMARY

Henley and Pollard⁽¹⁾ and Wroblewski and Ladue⁽²⁾ developed kinetic procedures for assaying alanine aminotransferase, formerly called glutamate pyruvate transaminase. The procedures were based on the oxidation of NADH by lactate dehydrogenase (LDH). This procedure incorporates significant modifications including optimization of substrate concentrations, and replacement of phosphate with Tris as buffer. These modifications meet the recommendations outlined by the International Federation of Clinical Chemistry (IFCC) for ALT measurement^(3,4). This procedure also uses stabilized NADH analog⁽⁵⁾.

Elevated activities of alanine aminotransferase result from a myocardial infarction, or hepatic disease. In cases of hepatic necrosis, elevation of ALT activities occur prior to the onset of clinical symptoms such as jaundice. To a lesser extent elevated activities are indicative of diseases of some internal organs⁽⁶⁾.

TEST PRINCIPLE

The ALT catalyzes the conversion of L-alanine and α -ketoglutarate to pyruvate and L-glutamate. In reaction II, LDH catalyzes the oxidation of the reduced cofactor to the cofactor. The rate of decrease in absorbance of the reaction mixture at 340 nm, due to the oxidation of the reduced cofactor, is directly proportional to the ALT activity.

REAGENTS

ALT L3K[®] Enzyme Reagent: A solution containing a buffer (pH 7.5 at 25°C), 500 mmol/L L-alanine, ≥ 12000 U/L lactate dehydrogenase (mammalian), 15 mmol/L α -ketoglutarate, 0.25 mmol/L NADH analog, stabilizers, and a preservative.

WARNINGS & PRECAUTIONS FOR USE

Avoid ingestion.

S24/25: Avoid contact with skin and eyes.

See Material Safety Data Sheet for additional information.

REAGENT STABILITY AND STORAGE

Reagents are ready for use.

Supplied reagent stable at 2-8°C until expiry date. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum.

SAMPLE STORAGE

Samples may be stored at 20-25°C for 24 hours or 2-8°C for 7 days.⁽⁷⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽¹⁰⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

This procedure measures total ALT. Red blood cells contain high concentrations of ALT; therefore, hemolysis can elevate results.

Interferences from icterus, lipemia and ascorbic acid were evaluated for this ALT method on a Roche/Hitachi 717[®] analyzer using a significance criterion of >10% variance from control.

Activity of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant	
42 U/L	Bilirubin	32 mg/dL	547 μ mol/L
43 U/L	Ascorbic Acid	2820 μ g/dL	160 μ mol/L
37 U/L	Intralipid	600 mg/dL	1800 mg/dL (20 mmol/L) Simulated Triglycerides

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁸⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' ALT L3K[®] reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material (If applicable).
3. Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:30, and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S, please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

The frequency of calibration, if necessary, using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the alanine aminotransferase activity of each sample.

TEST LIMITATIONS

A sample with an alanine aminotransferase level exceeding the linearity limit should be diluted with 0.9% saline and re-assayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽⁹⁾

Newborn: 13-45 U/L
Adult: Male: 10-40 U/L
Female: 7-35 U/L

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi 911[®] analyzer unless otherwise stated.

RESULTS

Alanine Aminotransferase activity is reported as U/L.

REPORTABLE RANGE (CLSI EP6)⁽¹⁰⁾

The linearity of the procedure described is 800 U/L. The lower limit of detection of the procedure described is 5 U/L. This data results in a reportable range of 5-800 U/L.

PRECISION STUDIES (CLSI EP5)⁽¹⁰⁾

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Within run precision was established by assaying each sample 20 times.

Activity (U/L)	Total SD (U/L)	Total CV (%)	Activity (U/L)	Within Run SD (U/L)	Within Run CV (%)
40	2.5	6.2	44	1.8	4.0
139	4.5	3.2	138	2.0	1.5

ACCURACY (CLSI EP9)⁽¹⁰⁾

The performance of this method (y) on a Roche/Hitachi 911[®] was compared with the performance of another ALT method (x) on a Roche/Hitachi 717[®]. Forty patient serum samples ranging from 7-469 U/L gave a correlation coefficient of 0.9994. Linear regression analysis gave the following equation:

$$\text{This method} = 0.994 (\text{reference method}) + 0.8 \text{ U/L}$$

TRADEMARKS

L3K is the registered trademark of Sekisui. All other trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

SEKISUI
DIAGNOSTICS

The Americas
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

International
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

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Temperature limitation
Límites de temperatura

REFERENCES/ REFERENCIAS

1. Henley, K.S., Pollard, H.M. A New Method for the Determination of Glutamic Oxalacetic and Glutamic Pyruvic Transaminase in Plasma, *J. Lab. and Clin. Med.*, 46, 785 (1955).
2. Wroblewski, F., Ladue, J.S. Serum Glutamic Pyruvic Transaminase in Cardiac and Hepatic Disease, *Proc. Soc. Exp. Biol. Med.* 91, 569 (1956).
3. Bergmeyer, H.U., and Horder, M. IFCC Methods for the Measurement of Catalytic Concentration of Enzymes, Part 3 IFCC Method For Alanine Aminotransferase, *J. Clin. Chem. Clin. Biochem.* 18, 521 (1980).
4. Lines, J.G., Dennis, P.M., deCediel, N., Khayat, M.H., Krawczynski, M.J., (Editors), AST and ALT Methods Approved, IFCC News, No. 40 (1985/2), p. 10.
5. US Patent #5,801,006.
6. Burtis, C.A., Ashwood, E.R. (Ed.), *Tietz Textbook of Clinical Chemistry*, W.B. Saunders Co., Toronto, p. 790-791 (1994).
7. Kaplan, Lawrence A. and Pesce, Amadeo J., (Ed), *Clinical Chemistry: Theory, Analysis, Correlation*, 3rd ed. Mosby-Year Book Inc., St. Louis, Missouri (1996), p. 518.
8. Young, D.S. *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, Third Edition (1990).
9. Burtis, C.A., Ashwood, E.R. (Ed.), *Tietz Textbook of Clinical Chemistry*, W.B. Saunders Co., Toronto, p.1800 (1999).
10. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.

IN34325-9
July 4, 2011