## BILIRUBIN-SL (TOTAL) ASSAY

CATALOGUE NO.: 243-10 SIZE: R1: 1 x 100 mL, R2: 1 x 25 mL

### INTENDED USE

For the IN VITRO quantitative measurement of Bilirubin (Total) in serum.

### TEST SUMMARY

Bilirubin is a bile pigment normally found in serum as a result of red cell destruction. It is a product of hemoglobin breakdown by the reticuloendothelial system and exists in two forms. Unconjugated (indirect) bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and excreted.

The elevation of total serum bilirubin may occur due to hemolytic processes, liver disease, or a disorder of the biliary tract.

Traditional methods of measuring bilirubin are based on the reaction of bilirubin with a diazo reagent to form the colored compound: azo-bilirubin. The diazo reaction can be accelerated by the addition of various chemicals. For example, Malloy-Evelyn<sup>(1)</sup> used ethanol, Jendrassik-Grof<sup>(2)</sup> used caffeine, and Walters-Gerarde<sup>(3)</sup> used DMSO. Modifications of these methods included the addition of surfactants as solubilizing agents.<sup>(4)</sup>

In this method, a 2,4-dichlorophenyl diazonium salt is used as the diazo reagent and the reaction is facilitated by the use of a surfactant.

### TEST PRINCIPLE

Total Bilirubin + 2,4-Dichlorophenyl diazonium salt  $\xrightarrow{\text{Surfactant}}$  Azobilirubin

Bilirubin (conjugated and unconjugated) couples with the diazo reagent in the presence of a surfactant to form azobilirubin. The increase in absorbance at 546 nm is directly proportional to the total bilirubin concentration.

### REAGENTS

Total Bilirubin-SL Accelerator Reagent R1: A solution containing 40 mmol/L phosphate buffer and 0.25 mol/L NaCl.

Total Bilirubin-SL Diazo Reagent R2: A solution containing 1 mmol/L 2,4-Dichlorophenyl diazonium salt and 30 mmol/L HCl.

### WARNINGS & PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.  
See Material Safety Data Sheet for additional information.

### REAGENT, PREPARATION, STORAGE & STABILITY

Reagents are ready for use.

Supplied reagents are stable at 2-8°C until expiry date.

Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum.

### SAMPLE STORAGE

Specimens must be protected from light. Samples may be stored at 2-8°C for 3 days or at minus 70°C for 3 months in the dark.<sup>(5)</sup>

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(7)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from hemoglobin and lipemia were evaluated for this total bilirubin method on a Roche/Hitachi<sup>®</sup> analyzer using a significance criterion of > 10% variance from control. One concentration of total bilirubin was evaluated. Interference data was collected in serum.

Hemoglobin produces significant interference with this method; therefore, hemolysed samples should not be used.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conv. Units	SI Units			
0.56 mg/dL	9.5 µmol/L	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

A summary of the influence of drugs on assays may be found by consulting Young, D.S.<sup>(6)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Sekisui Diagnostics' Total Bilirubin-SL reagents.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material.
3. Quality Control materials.

#### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:32:8 and a wavelength reading of 546 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration, if necessary, using an automated system is dependant on the system and the parameters used.

#### QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

#### CALCULATIONS

The analyzer automatically calculates the total bilirubin concentration of each sample.

#### TEST LIMITATIONS

A sample with a total bilirubin concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

#### REFERENCE INTERVALS<sup>(5)</sup>

0.2-1.0 mg/dL (3.4-17.1 µmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

#### PERFORMANCE CHARACTERISTICS

Data presented was collected on an automated analyzer unless otherwise stated.

#### RESULTS

Total Bilirubin concentration is reported as mg/dL (µmol/L).

#### REPORTABLE RANGE (CLSI EP6)<sup>(7)</sup>

The linearity of the procedure described is 30.0 mg/dL. The lower limit of detection of the procedure described is 0.0 mg/dL. This data results in a reportable range of 0.0-30.0 mg/dL (0-513 µmol/L).

## PRECISION STUDIES (CLSI EP5)<sup>(7)</sup>

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV %	Within Run SD		Within Run CV %
mg/dL	µmol/L	mg/dL	µmol/L		mg/dL	µmol/L	
0.8	14	0.03	0.5	4.2	0.03	0.5	3.7
8.9	152	0.37	6.3	4.2	0.15	2.6	1.7

## ACCURACY (CLSI EP9)<sup>(7)</sup>

The performance of this method (y) was compared with the performance of a similar method (x) on a Roche/Hitachi® analyzer. Sixty-five patient serum samples ranging from 0.1 to 10.6 mg/dL (2 to 181 µmol/L) gave a correlation coefficient of 0.9985. Linear regression analysis gave the following equation:

$$\text{This method} = 0.97 (\text{reference method}) + 0.12 \text{ mg/dL (2 } \mu\text{mol/L)}.$$

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**SEKISUI**  
DIAGNOSTICS

### The Americas

Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265

Fax: 902-628-6504

Email: questions@sekisuidiagnostics.com

peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

### International

Sekisui Diagnostics (UK) Limited  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, UK

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