CARBON DIOXIDE L3K[®] ASSAY

CATALOGUE NUMBER:	299-30	SIZE:	6 x 30 mL
	299-17		5 x 100 mL
	299-50		4 x 500 mL
	299-55		1 x 500 mL
	299-80		1 x 1000 mL

INTENDED USE


For the IN VITRO quantitative measurement of carbon dioxide in serum and plasma.

TEST SUMMARY

Elevated blood CO₂ is almost synonymous with respiratory acidosis. The latter is restricted to clinical conditions with a primary increase in carbon dioxide in the inspired air or increased metabolic production of carbon dioxide.


Decreased blood CO₂ is almost synonymous with respiratory alkalosis. The latter is restricted to clinical conditions with a primary decrease in carbon dioxide which can result from increased pulmonary ventilation due to mechanical ventilation or stimulation of the respiratory center.⁽¹⁾

Classic techniques for the measurement of carbon dioxide (CO₂) involve the addition of acid to liberate the carbon dioxide and the measurement of carbon dioxide thus released by manometric, volumetric, or titrimetric techniques. These procedures are both time consuming and cumbersome. Sekisui Diagnostics' Carbon Dioxide L3K[®] assay is an enzymatic procedure, employing phosphoenolpyruvate carboxylase (PEPC)⁽²⁾ and a stabilized NADH analog,⁽⁵⁾ which is easy to use and applicable to routine laboratory instrumentation.

TEST PRINCIPLE

PEPC catalyses the first reaction which produces oxaloacetate. In the presence of MDH, the reduced cofactor is oxidized by oxaloacetate. The decrease in concentration of the reduced cofactor is monitored between 405 and 415 nm and is proportional to the total carbon dioxide concentration in the sample.

PEPC is specific for the bicarbonate ion (HCO₃⁻) and its action disturbs the following equilibrium which results in conversion of the CO₂ to HCO₃⁻.

**REAGENTS**

Carbon Dioxide Reagent: A solution containing buffer (pH 7.5 at 25°C), 12.5 mmol/L PEP, >400 U/L PEPC (microbial), >4100 U/L malate dehydrogenase (mammalian), 0.6 mmol/L NADH analog, activators, stabilizers, a surfactant, and a preservative.

WARNINGS & PRECAUTIONS FOR USE

S20: When using, do not eat or drink.
S24/25: Avoid contact with skin and eyes.
See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use.
Supplied reagent stable at 2-8°C until expiry date.
Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum or lithium heparinized plasma. The specimen should be promptly separated from the clot or cells and stored tightly sealed to prevent the loss of carbon dioxide.

SAMPLE STORAGE

Samples should be stored at 2-8°C.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁴⁾

Cross Contamination studies have not been performed on automated instruments. Certain reagent / instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia and hemolysis were evaluated for this carbon dioxide method on a Roche/Hitachi[®] 717 analyzer using a significance criterion of >10% variance from control. Data for plasma is similar to that presented for serum.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
26.3 mmol/L	26.3 mEq/L	Hemoglobin	600 mg/dL	93 µmol/L
27.1 mmol/L	27.1 mEq/L	Bilirubin	40 mg/dL	684 µmol/L
24.6 mmol/L	24.6 mEq/L	Intralipid	1000 mg/dL	3000 mg/dL (34 mmol/L) Simulated Triglycerides

The information presented above is based on results from Sekisui Diagnostic's studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁵⁾

ANALYTICAL PROCEDURE**MATERIALS PROVIDED**

Sekisui Diagnostics' Carbon Dioxide L3K[®] Reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a fixed rate test mode, with a sample to reagent ratio of 1:100, and a wavelength reading of 405 nm or 415 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the carbon dioxide concentration of each sample.

TEST LIMITATIONS

A sample with a carbon dioxide concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS

22-29 mmol/L (22-29 mEq/L)⁽⁶⁾

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on Roche/Hitachi[®] 717 analyzer unless otherwise stated.

RESULTS

Carbon dioxide concentration is reported as mmol/L (mEq/L).

REPORTABLE RANGE (CLSI EP6)⁽⁴⁾

The linearity of the procedure described is 50.0 mmol/L (mEq/L). The lower limit of detection of the procedure described is 2.9 mmol/L (mEq/L). This data results in a reportable range of 2.9 to 50.0 mmol/L (mEq/L).

PRECISION STUDIES (CLSI EP5)⁽⁴⁾

Total precision data was collected on two concentrations of control sera in triplicate over five runs, for a total of 15 sample values.

Concentration	Total SD	Total CV%	Concentration	Within Run SD	Within Run CV%
mmol/L (mEq/L)	mmol/L (mEq/L)		mmol/L (mEq/L)	mmol/L (mEq/L)	
14	0.24	1.7	14	0.18	1.3
25	0.27	1.1	24	0.40	1.7

Within run precision data was collected on two concentrations of control sera, each run 20 times in a single assay.

ACCURACY (CLSI EP9)⁽⁴⁾

The performance of this method (y) was compared with the performance of a similar carbon dioxide method (x) on a Cobas Mira[®] analyzer.

Forty-five serum samples ranging from 15.4 to 43.7 mmol/L (mEq/L) were tested and gave a correlation coefficient of 0.9861. Linear regression analysis gave the following equation:

$$\text{This method} = 1.00(\text{reference method}) - 0.42 \text{ mmol/L (mEq/L)}.$$

The performance of this method (y) with plasma samples (y) was compared to this method with serum samples (x) on an Advia® 1650 analyzer.

Twenty-five patient serum and plasma samples ranging from 17.1 to 25.3 mmol/L (mEq/L) gave correlation coefficient of 0.9636. Linear regression analysis gave the following equation:

$$\text{This method (plasma)} = 0.950[\text{this method(serum)}] + 1.35 \text{ mmol/L (mEq/L)}.$$

TRADEMARKS

L3K is a registered trademark of Sekisui. All other trademarks, brands, product names are the property of their respective companies.

Manufactured by:

The Americas

Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
peidiagnosticttechnical@sekisuidiagnostics.com
www.sekisuidiagnostics.com

International

Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK
Email: info@sekisuidiagnostics.com

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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REFERENCES/RÉFÉRENCES/REFERENCIAS/ RIFERIMENTI/ LITERATURNACHWEIS/ LITERATUUR

1. Tietz, N.W. (Ed.), *Fundamentals of Clinical Chemistry*, W.B. Saunders Co., Toronto, 636-638, 937 (1970).
2. Norris, K.A., Atkinson, A.R., Smith, W.G., Clin. Chem. 21 (1975).
3. US Patent #5,801,006
4. *CLSI Guidelines and Standards*, Clinical and Laboratory Standards Institute, Wayne, PA.
5. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Third Edition, Washington (1990).
6. Tietz, N.W. (Editor), *Textbook of Clinical Chemistry*, W.B. Saunders Co., Philadelphia (1986).

PATENT NUMBER
US Patent No. 5,801,006

Authorized Representative/Nom et adresse du représentant autorisé/Representante autorizado/Nome e indirizzo del rappresentante autorizzato/Name and Adresse des zugelassenen Vertreters/Naam en adres van gevlmachtigde vertegenwoordiger:
Sekisui Diagnostics (UK) Limited

50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF
United Kingdom
Tel (+44)(0)1732-220022
Fax (+44)(0)1732-220024

IN29930-15
June 14, 2011