WIDE RANGE CRP ASSAY

CATALOGUE NUMBER: 250-20 **SIZE:** R1: 1 x 20 mL, R2: 1 x 20 mL
250-25 R1: 1 x 50 mL, R2: 1 x 50 mL

INTENDED USE

For the IN VITRO quantitative measurement of C-Reactive Protein (CRP) in serum.

TEST SUMMARY

CRP is one of the cytokine-induced “acute phase” proteins whose blood concentration rises during a general and unspecific response to infections and non-infectious inflammatory processes. Intra-individual variations of CRP are significant and should be taken into account when interpreting values. They can range from 30-60%.⁽¹⁾ Because the increase is nonspecific, it cannot be interpreted without a complete history.⁽²⁾

Past methods for the measurement of CRP included semi-quantitative analysis by precipitation with C-polysaccharide in the presence of calcium ions.

Immunological methods have since been developed to allow quantitative analysis of CRP. This method, using antihuman CRP mouse antibody coated latex to bind CRP, allows for accurate measurement over a wide range of CRP concentrations. The resulting agglutination from binding is measured spectrophotometrically and is proportional to the concentration of CRP in the sample.

TEST PRINCIPLE

CRP + antihuman CRP mouse monoclonal antibody-coated latex → immune complexes

CRP reacts with antihuman CRP mouse antibody-coated latex, and agglutination occurs. CRP concentration is determined by measurement of the change of absorbance that results from the agglutination reaction.

REAGENTS

Wide Range CRP Buffer Solution (R1): A solution containing a buffer at pH 8.5.

Wide Range CRP Latex Reagent (R2): A solution containing antihuman CRP mouse monoclonal antibody-coated latex (2 mg/dL).

WARNINGS & PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.
See Material Safety Data Sheet for additional information.
This product has not been tested for HAMA (human anti-mouse antibodies).

REAGENT PREPARATION, STORAGE, AND STABILITY

Reagents are ready for use.
Supplied reagents are stable at 2-8°C until the expiry date stated on the vial label.
Stability claims are based on real time studies.

REAGENT DETERIORATION

The buffer solution should be clear. Turbidity would indicate deterioration.

The latex reagent is cloudy. Values of control material outside acceptable range may be an indication of instability.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum.

SAMPLE STORAGE

Samples should be analyzed fresh, stored at 4°C for < 72 hours, -20°C for six months or -70°C indefinitely.⁽⁶⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁵⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia, and hemolysis and ascorbic acid were evaluated for this method on a Roche/Hitachi® 911 analyzer using a significance criterion of >10% variance from control.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
mg/dL	mg/L			
0.3	3	Hemoglobin	1000 mg/dL	155 µmol/L
0.3	3	Bilirubin	50 mg/dL	855 µmol/L
0.3	3	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides
0.3	3	Ascorbic Acid	100 mg/dL	5.7 mmol/L

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽³⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' Wide Range CRP reagents.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbances at appropriate wavelength as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:50:50, and wavelength readings of (primary/secondary) 570 nm / 800 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. Sekisui's Wide Range CRP Multicalibrator Set, Cat. No. SE-250, is recommended. The frequency of calibration, if necessary, is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the C-Reactive Protein concentration of each sample.

TEST LIMITATIONS

A sample with a C-Reactive Protein concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed, incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽⁴⁾

Less than 0.3 mg/dL (3 mg/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 911 analyzer unless otherwise stated.

RESULTS

C-Reactive Protein concentration is reported as mg/dL (mg/L).

REPORTABLE RANGE (CLSI EP6)⁽⁵⁾

The reportable range is dependent upon the calibrators used to calculate the non-linear calibration curve. The lower limit of detection for the procedure described is 0.0 mg/dL (0 mg/L). The lower limit of quantitation for the procedure described is 0.0 mg/dL (0 mg/L). The highest calibrator for the procedure described is 36.0 mg/dL (360 mg/L), resulting in a reportable range of 0.0 to 36.0 mg/dL (0 to 360 mg/L). CRP concentrations up to at least 100 mg/dL (1000 mg/L) do not demonstrate prozone effects.

PRECISION STUDIES (CLSI EP5)⁽⁵⁾

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV %	Concentration		Within Run SD		Within Run CV %
mg/dL	mg/L	mg/dL	mg/L		mg/dL	mg/L	mg/dL	mg/L	
0.2	2	0.00	0.10	3.7	0.2	2	0.00	0.00	1.9
1.6	16	0.04	0.38	2.3	1.6	16	0.02	0.23	1.5

Within run precision data was collected on two concentrations of control sera each run 20 times in a single assay.

ACCURACY (CLSI EP9)⁽⁵⁾

The performance of this method (y) was compared with the performance of a similar method (x) on a Roche/Hitachi® 911. Forty patient serum samples ranging from 0.1 to 27.4 mg/dL (1 to 274 mg/L) were tested and gave a correlation coefficient of 0.9999. Linear regression analysis gave the following equation:

$$\text{This method} = 0.975 (\text{reference method}) + 0.017 \text{ mg/dL (0.17 mg/L)}$$

The performance of this method (y) was compared with the performance of a similar method (x) on a Roche/Hitachi® 911. Forty patient serum samples ranging from 0.1 to 0.9 mg/dL (1 to 9 mg/L) were tested and gave a correlation coefficient of 0.9979. Linear regression analysis gave the following equation:

$$\text{This method} = 0.992 (\text{reference method}) - 0.0027 \text{ mg/dL (0.027 mg/L)}$$

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

The Americas
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

International
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
peidiagnosticttechnical@sekisuidiagnostics.com

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Temperature limitation
Límites de temperatura

REFERENCES/ REFERENCIAS

1. Callaghan, J.V., Guidance for Industry and FDA Staff - Review Criteria for Assessment of C-Reactive Protein (CRP), High Sensitivity C-Reactive Protein (hsCRP) and Cardiac C-Reactive Protein (cCRP) Assays, Issue Date September 22, 2005.
2. Burtis, C.A., Ashwood, E.R. (Ed.), Tietz Textbook of Clinical Chemistry, W.B. Saunders Co., Toronto, Third Edition p. 493.
3. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, Third Edition.
4. Liuzzo, G., et al., N. Eng. J. Med., 331:417-424 1994.
5. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
6. Wu, Alan H.B., et al., Tietz Clinical Guide to Laboratory Tests, 4th Edition, p. 190.

IN25020-6
August 16, 2011