## TOTAL PROTEIN (MICRO) ASSAY

**CATALOGUE NUMBER:** 450-50 **SIZE:** 2 x 250 mL + 15 mL Calibrator

### INTENDED USE

For the IN VITRO quantitative measurement of total protein concentration in urine and cerebrospinal fluid (CSF).

### TEST SUMMARY

The presence of protein in the urine is useful in the evaluation of renal disorders. The degree of proteinuria varies with the type and severity of the disease and measuring the amount of protein excreted aids in the diagnosis and directs the therapy. Excess protein in CSF is indicative of infection and its measurement is useful in diagnosis.

This method for the determination of urine proteins is based on the procedure developed by Fujita<sup>(1)</sup> and Watanabe<sup>(2)</sup>. It is a sensitive dye-binding colorimetric method. Pyrogallol red is combined with molybdenum acid. At a low pH, the dye is red, the color changes to blue under test conditions when complexed with proteins. This method reduces cuvette staining and is easily adapted to automated analyzers.

### TEST PRINCIPLE

Protein + Pyrogallol → Protein Pyrogallol  
Red-Molybdenum Complex Red-Molybdenum Complex

The increase in absorbance at 600 nm due to the formation of the colored complex is directly proportional to the concentration of protein in the reaction.

### REAGENTS

**Total Protein (Micro) Reagent:** An aqueous solution containing a buffer, 0.067 mmol/L pyrogallol red, at least 1 mmol/L molybdenum, a stabilizer, a surfactant, and a preservative.

**Total Protein (Micro) Calibrator:** An aqueous solution containing 100 mg/dL (1 g/L) human serum albumin and a preservative.

### WARNINGS & PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.

The reagent provided is for IN VITRO diagnostic use. Do not pipette the reagent by mouth and avoid contact with eyes.

Any donor unit which may be used in the preparation of this product has been tested and found non-reactive for Hepatitis B surface antigen (HBsAg) and HIV antibody. No known test method can offer complete assurance that products derived from human blood will not transmit infectious agents. Therefore, all human based blood products and patients samples are considered potentially hazardous and should be handled in accordance with good laboratory practices.

See Material Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use.

Supplied reagents are stable at 2-8°C until expiry date.

Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed urine or cerebrospinal fluid (CSF). Filter or centrifuge turbid samples to remove cells. Urine samples may be random or 24 hour collections. Do not use preservatives.

### SAMPLE STORAGE

Samples of urine and CSF should be assayed promptly. Samples may be stored frozen for several weeks if further testing is necessary.

### ANALYTICAL SPECIFICITY

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

The presence of hemoglobin in urine and CSF will falsely elevate results. Urine colored as a result of drugs and xanthochromic CSF require the use of a sample blank. A blank may be prepared using 20 µL of sample and 2.0 mL of deionized water and should be read at 600 nm.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(3)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Sekisui Diagnostics' total protein (micro) reagent and calibration material.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
- 2) Quality Control materials.

#### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed using an endpoint test mode, with a sample to reagent ratio of 1:100, and a wavelength reading of 600 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependant on the system and parameters used.

#### QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state, and federal regulations. The results should fall within the acceptable range as established by the laboratory.

#### CALCULATIONS

The analyzer automatically calculates the total protein (micro) concentration of each sample.

#### TEST LIMITATIONS

A sample with a total protein (micro) concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

#### REFERENCE INTERVALS

#### EXPECTED VALUES

Urine:	28-141 mg/day <sup>(2)</sup>
CSF Adults:	15-40 mg/dL (0.15-0.40 g/L) <sup>(4)</sup>
CSF Newborn, < 1 month:	20-120 mg/dL (0.20-1.20 g/L) <sup>(4)</sup>

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

#### PERFORMANCE CHARACTERISTICS

#### RESULTS

Total protein (micro) concentration is reported as mg/dL (g/L).

#### REPORTABLE RANGE

The linearity of the procedure described is 200 mg/dL (2 g/L). Linearity will depend on the sample/ reagent ratio used.

## PRECISION STUDIES

Within day precision data was collected on two concentrations of urine controls. Nineteen replicates of level 1 and thirty replicates of level 2 were assayed.

Total Protein (Micro)	Concentration		Standard Deviation		Coefficient of Variation %
	mg/dL	g/L	mg/dL	g/L	
Level 1	16.0	0.16	0.90	0.009	5.5
Level 2	40.0	0.40	1.40	0.014	3.6

Day to day precision data was collected on two concentrations of CSF controls assayed ten times over ten days.

Total Protein (Micro)	Concentration		Standard Deviation		Coefficient of Variation %
	mg/dL	g/L	mg/dL	g/L	
Level 1	34.0	0.34	2.30	0.023	6.9
Level 2	65.0	0.65	3.20	0.032	4.9

Run to run precision data was collected on two concentrations of urine controls assayed 15 times over five runs.

Total Protein (Micro)	Concentration		Standard Deviation		Coefficient of Variation %
	mg/dL	g/L	mg/dL	g/L	
Level 1	17.0	0.17	2.30	0.023	13.5
Level 2	56.0	0.56	3.30	0.033	5.9

Run to run precision data was collected on two concentrations of CSF assayed 14 times over seven runs.

Total Protein (Micro)	Concentration		Standard Deviation		Coefficient of Variation %
	mg/dL	g/L	mg/dL	g/L	
Level 1	40.0	0.40	2.40	0.024	6.2
Level 2	61.0	0.61	1.40	0.014	2.4

## ACCURACY

The performance of this method (y) was compared with the performance of a similar method (x). Forty-three urine samples ranging from 21.0-138.0 mg/dL (0.21-1.38 g/L) gave a correlation coefficient of 0.9826. Linear regression analysis gave the following equation:

$$\text{This method} = 1.12 (\text{reference method}) - 3.0 \text{ mg/dL} (0.03 \text{ g/L}).$$

The performance of this method (y) was compared with the performance of a similar method (x). Forty CSF samples ranging from 18.0-168.0 mg/dL (0.18-1.68 g/L) gave a correlation coefficient of 0.9879. Linear regression analysis gave the following equation:

$$\text{This method} = 1.078 (\text{reference method}) - 1.9 \text{ mg/dL} (0.019 \text{ g/L}).$$

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**SEKISUI**  
DIAGNOSTICS

**The Americas**  
Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265  
Fax: 902-628-6504

Email: questions@sekisuidiagnostics.com  
peidiagnostictchnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

**International**  
Sekisui Diagnostics (UK) Limited  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, UK

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Temperature limitation  
Límites de temperatura

## REFERENCES/ REFERENCIAS

1. Fujita, Y., Mori, I., and Kitano, S., *Bunseki, Kagaku*, 32, 379 (1983).
2. Watanabe, N., Kamel, S., Ohkubo, A., Tokudo, K., *Urinary Protein as Measured with a Pyrogallol Red-Molybdate Complex, manually and in a Hitachi 726 automated analyzer*, *Clin. Chem.* 32, 8, 1551 (1986).
3. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, 1990.
4. Tietz, N.W., *Textbook of Clinical Chemistry*, W.B. Saunders Co. 1986, p. 1844.

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Sekisui Diagnostics (UK) Limited  
50 Gibson Drive  
Kings Hill, West Malling  
Kent, ME19 4AF  
United Kingdom/ Reino Unido  
Tel (+44)(0)1732-220022  
Fax (+44)(0)1732-220024

IN45050-15  
August 16, 2011