osom[®]

QC Inside

Mono Test

EN	1
DK	6
FR	11
DE	16
IT	21
NO	26
ES	31
SV	36

**DEFINITIONS OF SYMBOLS / SYMBOLFORKLARING / DÉFINITIONS DES SYMBOLES /
DEFINITIONEN DER SYMBOLE / DEFINIZIONE DEI SIMBOLI / SYMBOLER /
DEFINICIONES DE SÍMBOLOS / SYMBOLDEFINITIONER**

..... 41

**REFERENCES / REFERENCER / RÉFÉRENCES / LITERATUR / BIBLIOGRAFIA /
REFERANSER / BIBLIOGRAFÍA / REFERENSER**

..... 42

EN

OSOM[®] Mono Test

Catalog number 145E

CLIA Complexity: Non-Waived for Serum and Plasma

CLIA Complexity: Waived for Whole Blood

FOR LABORATORY AND PROFESSIONAL IN VITRO DIAGNOSTIC USE ONLY.

INTENDED USE

The OSOM[®] Mono Test is intended for the qualitative detection of infectious mononucleosis heterophile antibodies in serum, plasma or whole blood as an aid in the diagnosis of infectious mononucleosis.

SUMMARY AND EXPLANATION OF TEST

The diagnosis of infectious mononucleosis (IM) is suggested on the basis of the clinical symptoms of fever, sore throat and swollen lymph glands. The highest incidence of symptomatic IM occurs during late adolescence (15-24 years of age). Infectious mononucleosis is caused by the Epstein-Barr Virus (EBV)^{1,2}. The laboratory diagnosis of IM is based on the detection of IM heterophile antibodies. These heterophile antibodies are directed against antigens found in bovine, sheep and horse erythrocytes. The OSOM Mono Test utilizes an extract of bovine erythrocytes to give the required sensitivity and specificity.

PRINCIPLES OF TEST

The OSOM Mono Test uses color immunochromatographic dipstick technology with bovine erythrocyte extract coated on the membrane. In the test procedure, serum, plasma or whole blood is mixed with the Diluent. Then the Test Stick is placed in the mixture and the mixture migrates along the membrane. If the specific IM heterophile antibody is present in the sample, it will form a complex with the bovine erythrocyte extract conjugated color particles. The complex will then be bound by bovine erythrocyte extract immobilized on the membrane and a visible blue Test Line will appear to indicate a positive result.

KIT CONTENTS AND STORAGE

25 Test Sticks in a container

25 Test Tubes

25 Transfer Pipettes

25 Capillary Tubes with 1 Capillary Bulb

1 Diluent (contains buffer with 0.2% sodium azide)

1 Mono Positive Control (contains rabbit anti-beef stroma in tris buffer with 0.2% sodium azide and 0.05% gentamycin sulfate preservatives)

1 Mono Negative Control (contains goat albumin in tris buffer with 0.2% sodium azide)

1 Work Station

1 Directional Insert

2 Additional test sticks have been included in the kit for external QC testing

Note: Extra components (tubes, pipettes, capillary tubes, capillary bulb) have been provided for your convenience.

Store the Test Sticks and reagents tightly capped at 15°-30°C (59°-86°F).

Do not use the Test Sticks or reagents after their expiration dates.

MATERIALS REQUIRED BUT NOT PROVIDED

Specimen collection containers.

A timer or watch.

PRECAUTIONS

- For in-vitro diagnostic use only.
- Follow your laboratory safety guidelines in the collection, handling, storage and disposal of patient specimens and all items exposed to patient specimens.
- The Diluent and Controls contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azide. Large quantities of water must be used to flush discarded Diluent or Controls down a sink.
- The Capillary Bulb contains dry natural rubber.
- Do not interchange or mix components from different kit lots.

The Diluent, Positive Control and Negative Control contain 0.2% Sodium Azide.

In the EU this is classified as Harmful with the following Risk and Safety phrases:

R22 Harmful if swallowed

R32 Contact with acids liberates very toxic gas

S35 This material and its container must be disposed of in a safe way.

SPECIMEN COLLECTION AND PREPARATION

Serum, Plasma or Whole Blood Sample

Obtain specimens by acceptable medical technique. Collect whole blood samples using a tube containing EDTA or heparin as an anticoagulant. Other anticoagulants have not been tested. Serum and plasma specimens may be refrigerated (2°-8°C; 36°-46°F) and tested within 48 hours; serum and plasma specimens held for longer times should be frozen (below -10°C; 14°F) and tested within 3 months. Test whole blood specimens within 24 hours. Specimens must be at room temperature (15°-30°C; 59°-86°F) when tested.

Fingertip Whole Blood

Hold the capillary tube horizontally while collecting the sample. Touch the end of the capillary tube to the drop of blood on the patient's finger. Fill the capillary tube completely. Place the small end of the black bulb onto the capillary tube. Place your fingertip over the opening in the bulb. Squeeze the bulb to dispense the whole blood sample into the test tube.

QUALITY CONTROL

External Quality Control

For external QC testing, use the controls provided in the kit. Add one free falling drop of control to the Test Tube and then proceed in the same manner as with a patient sample. Quality Control requirements should be established in accordance with local, state and federal regulations or accreditation requirements. Minimally Sekisui Diagnostics recommends that positive and negative external controls be run with each new lot and with each new untrained operator. Some commercial controls may contain interfering additives. The use of these controls is not recommended.

Internal Quality Controls

The OSOM Mono Test provides two levels of internal procedural controls with each test procedure.

- The red Control Line is an internal positive control. The Test Stick must absorb the proper amount of sample and be working properly for the red Control Line to appear.
- A clear background is an internal negative control. If the test has been performed correctly and the Test Stick is working properly, the background will clear to give a discernible result.

If the red Control Line does not appear, the test is invalid. If the background does not clear and interferes with the test result, the test may be invalid.


LIMITATIONS

- As with all diagnostic assays, the results obtained by this test yield data that must be used as an adjunct to other information available to the physician.
- The OSOM Mono Test is a qualitative test for the detection of IM heterophile antibody.
- A negative result may be obtained from patients at the onset of the disease due to heterophile antibody levels below the sensitivity of this test kit. If symptoms persist or intensify, the test should be repeated.
- Some segments of the population with acute IM are heterophile antibody negative¹.

TEST PROCEDURE

- Addition of Specimen:

- Slowly add 1 drop of Diluent to the bottom of the Test Tube.

- Mix.

- Remove the Test Stick(s) from the container. Re-cap the container immediately.

- Place the Absorbent End of the Test Stick into the treated sample. Leave the Test Stick in the Test Tube.

- Read results at 5 minutes. Positive results may be read as soon as the red Control Line appears.


INTERPRETATION OF TEST RESULTS

Notes


A blue or red line which appears uneven in color density is considered a valid result.

Positive


A blue Test Line and a red Control Line is a positive result for the detection of infectious mononucleosis heterophile antibody. Note that the blue line can be any shade of color and can be lighter or darker than the line in the picture.

Negative

A red Control Line but no blue Test Line is a negative result. No infectious mononucleosis heterophile antibody has been detected.

Invalid

If after 5 minutes, no red Control Line appears or background color makes reading the red Control Line impossible, the result is invalid. If this occurs, repeat the test on a new Test Stick.

EXPECTED VALUES

A heterophile antibody response is observed in approximately 80-90% of adults and children with EBV-caused IM. This percentage drops to approximately 50% for children under four years of age¹.

While the incidence of IM reflects wide seasonal, ethnic and geographical variation, a large epidemiological study noted that the highest incidence of symptomatic IM occurred during late adolescence (15-24 years of age)².

PERFORMANCE CHARACTERISTICS

A total of 439 specimens (183 serum, 176 plasma and 80 whole blood) were evaluated by two clinical labs in a clinical study. Test results of the OSOM Mono Test were compared to results obtained with a commercially available latex particle agglutination test for the qualitative determination of infectious mononucleosis heterophile antibodies. Discrepancies between the results given by the OSOM Mono Test and the latex particle agglutination test were resolved by Epstein-Barr Virus (EBV) specific serological assays. In these assays, the specific antibodies to the EBV capsid antigen (IgM) and EBV nuclear antigen-1 (IgM and IgG) were determined.

Serum Specimens:		Comparative Test	
		+	-
	OSOM Mono Test	+ 74	8*
	-	0	101

*6 out of 8 tested positive by EBV testing

Plasma Specimens:		Comparative Test	
		+	-
	OSOM Mono Test	+ 67	15*
	-	0	94

*8 out of 15 tested positive by EBV testing

Whole Blood Specimens:		Comparative Test	
		+	-
	OSOM Mono Test	+ 30	3*
	-	0	47

*1 out of 3 tested positive by EBV testing

All Specimens:		Comparative Test	
		+	-
	OSOM Mono Test	+ 171	26*
	-	0	242

*15 out of 26 tested positive by EBV testing

When compared to a commercially available latex particle agglutination test for infectious mononucleosis heterophile antibodies, the OSOM Mono Test showed a sensitivity of 100% and a specificity of 90.3%. The overall agreement was 94.1%.

Fifteen of the twenty-six discrepant samples were determined to be recent or acute EBV infections by EBV serological testing, in which case the sample was considered positive. Including the samples confirmed positive by EBV serological testing, the overall clinical study specificity of the OSOM Mono Test is 95.9% and the overall sensitivity is 100%.

POL Studies

An evaluation of the OSOM Mono Test was conducted at three physicians' offices or clinical laboratories where testing was performed by personnel with diverse educational backgrounds. Each site tested the randomly coded panel consisting of negative (5), low positive (3) and moderate positive (4) specimens for three days. The results obtained had 99.1% agreement (107/108) with the expected results.

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IVD

In vitro diagnostic medical device

Medicinsk udstyr til in vitro-diagnostik
Dispositif médical de diagnostic in vitro
In-Vitro-Diagnostikum
Dispositivo medico-diagnostico in vitro
Til in vitro-diagnostisk bruk
Producto sanitario para diagnóstico in vitro
Medicintekniske produkter for in vitro diagnostik

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EC REP

Authorized representative in the EC

Repræsentant i det Europæiske Fællesskab
Mandataire dans la Communauté européenne
Bevollmächtigter in der Europäischen Gemeinschaft
Mandatario nella Comunità Europea
Autorisert representant
Representante autorizado en la Comunidad Europea
Auktoriserad representant i Europeiska Gemenskapen

CONTROL +

Positive control

Positiv kontrol
Contrôle positif
Positiv kontrolle
Controllo positivo
Positiv kontroll
Control positivo
Positiv kontroll

Consult instructions for use

Se brugsanvisning
Consulter les instructions d'utilisation
Gebrauchsanweisung beachten
Consultare le istruzioni per l'uso
Se bruksanvisningen
Consultar las instrucciones de uso
Se handhavandebeskrivningen

CONTROL -

Negative control

Negativ kontrol
Contrôle négatif
Negative kontrolle
Controllo negativo
Negativ kontroll
Control negativo
Negativ kontroll

Caution, consult accompanying documents

Forsigtig se brugsanvisning
Attention voir notice d'instructions
Achtung, Begleitdokumente beachten
Attenzione, vedere le istruzioni per l'uso
Forsiktig, se vedlagte dokumenter
Atención, ver instrucciones de uso
Försiktighet, se handhavandebeskrivningen

Temperature limitation

Temperaturbegrensning
Limites de température
Temperaturbegrenzung
Limiti di temperatura
Temperaturbegrensning
Limite de temperatura
Temperaturbegrensning

Contains sufficient for "n" tests

Indeholder tilstrækkeligt til "n" test
Contenu suffisant pour <n> tests
Inhalt ausreichend für <n> Prüfungen
Contenuto sufficiente per "n" saggi
Antall tester
Contenido suficiente para <n> ensayos
Räcker till "n" antal tester

**REFERENCES / REFERENCER / RÉFÉRENCES / LITERATUR / BIBLIOGRAFIA
REFERANSER / BIBLIOGRAFÍA / REFERENSER**

1. Lennette, E.T., Epstein-Barr Virus, in Manual of Clinical Microbiology. Balows, A., Hausler, W.J. Jr., Hermann, K.L., Isenberg, H.D., Shadomy, H.J., Editors, 5th Edition, American Society for Microbiology, Washington D.C., pp847-852, 1991.
2. Heath, C.W. Jr., Brodsky, A.L., Potolsky, A.I., Infectious Mononucleosis in a General Population. Am. J. Epidemiol., 95:46, 1972.

OSOM® and the QC Inside logo are registered U.S. trademarks of Sekisui Diagnostics, LLC er registrerede amerikanske varemærker tilhørende Sekisui Diagnostics, LLC sont des marques déposées américaines de Sekisui Diagnostics, LLC sind eingetragene US-Marken von Sekisui Diagnostics, LLC sono marchi registrati degli Stati Uniti di Sekisui Diagnostics, LLC er registrerte varemerker i USA for Sekisui Diagnostics, LLC son marcas registradas de EE.UU. Sekisui Diagnostics, LLC är registrerade varumärken i USA av Sekisui Diagnostics, LLC

Licensed under U.S. Patent Nos. 5,714,389; 5,989,921 and 6,485,982 and related non-U.S. patents and patent applications.

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Sekisui Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF, UK

SEKISUI
DIAGNOSTICS

Framingham, MA 01701 USA
www.sekisuidiagnostics.com

OSOM[®] Mono Test

1

For serum, plasma, or whole blood Samples in tubes

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