SALICYLATE-SL ASSAY

CATALOGUE NUMBER: 511-20 **SIZE:** 2 x 10 mL, 2 x 10 mL
511-40 4 x 10 mL, 4 x 10 mL

Note: Changes are highlighted.

INTENDED USE

For the IN VITRO quantitative measurement of salicylate concentration in serum and plasma.

TEST SUMMARY


Salicylate is a common drug used for its analgesic and anti-inflammatory properties. Its accessibility leads to its implication in a large number of accidental ingestions by children and it is a common choice among adults and adolescents for attempted suicidal poisoning.⁽¹⁾

Salicylate overdose results in disturbances of the central nervous system and the gastrointestinal tract as well as encephalopathy and renal failure. It represents an acute medical emergency and rapid quantitation of the drug is necessary for effective patient management.⁽²⁾

Salicylate has been traditionally measured by the "Trinder Reaction" which is based on the interaction between salicylate and ferric ions. This test however is not specific and requires extraction or centrifugation which inhibit the automation of the test.

This enzymatic Salicylate Assay provides a rapid, specific and simplified method for salicylate measurement. It is based on the action of salicylate hydroxylase on salicylate and NADH which results in a decrease in absorbance proportional to the amount of salicylate present.⁽³⁾ The test can be adapted to automated instruments resulting in rapid, accurate results required by physicians.

TEST PRINCIPLE

Salicylate hydroxylase catalyzes the conversion of salicylate and NADH to catechol and NAD in the presence of oxygen. The resulting decrease in absorbance at 340 nm, due to the conversion of NADH to NAD, is directly proportional to the concentration of salicylate in the sample.

REAGENTS

Salicylate-SL Enzyme Reagent: A solution containing a buffer (pH 7.6 at 25°C), at least 770 U/L salicylate hydroxylase (microbial), stabilizers and a preservative.

Salicylate-SL NADH Reagent: A solution containing a buffer (pH 10.4 at 25°C) and 1.2 mmol/L NADH.

Salicylate Calibrator: 1 x 5 mL of a solution containing 20.7 mg/dL (1.50 mmol/L) salicylate, and a preservative.

WARNINGS & PRECAUTIONS FOR USE

S20: When using, do not eat or drink.
S24/25: Avoid contact with skin and eyes.
See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use. Supplied reagents are stable at 2-8°C until expiry date. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum or lithium heparinized plasma.

SAMPLE STORAGE

Specimens for analysis should be assayed promptly.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁴⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia, and hemolysis were evaluated for this salicylate method on a Roche/Hitachi® 717 analyzer using a significance criterion of >10% variance from control. Interference data was collected in serum. Plasma data is expected to be similar.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
35.6 mg/dL	2.58 mmol/L	Hemoglobin	400 mg/dL	62 µmol/L
36.9 mg/dL	2.67 mmol/L	Bilirubin	40 mg/dL	684 µmol/L
36.5 mg/dL	2.64 mmol/L	Intralipid	800 mg/dL	2400 mg/dL (27.1 mmol/L) Simulated Triglycerides

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁵⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' Salicylate reagents and Salicylate calibrator.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
- Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:16:33(Dil):35 and wavelength readings of (primary/secondary) 340/376 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

A salicylate calibrator is included with the reagents and should be used as directed to calibrate the procedure. The frequency of calibration on automated systems is dependent on the system and the parameters used.

QUALITY CONTROL

Appropriate concentrations of quality control material should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the salicylate concentration of each sample.

TEST LIMITATIONS

A sample with a salicylate concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽⁶⁾

Toxic concentration: greater than 30.0 mg/dL (2.17 mmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 717 analyzer unless otherwise stated.

RESULTS

Salicylate concentration is reported as mg/dL (mmol/L).

REPORTABLE RANGE (CLSI EP6)⁽⁴⁾

The linearity of the procedure described is 100.0 mg/dL (7.24 mmol/L). The limit of quantitation of the procedure described is 5.3 mg/dL (0.38 mmol/L). This data results in a reportable range of 5.3-100.0 mg/dL (0.38-7.24 mmol/L).

PRECISION STUDIES (CLSI EP5)⁽⁴⁾

Total precision data was collected on four concentrations of control sera using a single lot of reagent in 40 runs conducted over 20 days. Within run precision data was collected by assaying twenty samples of four concentrations of control sera in one run using one lot of reagent.

Concentration		Total SD		Total CV%	Concentration		Within Run SD		Within Run CV%
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	mg/dL	mmol/L	
9.2	0.67	0.66	0.048	7.2	8.9	0.64	0.43	0.031	4.8
16.6	1.20	0.90	0.065	5.4	16.3	1.18	0.60	0.043	3.7
34.9	2.53	1.18	0.085	3.4	33.5	2.43	0.52	0.038	1.5
49.5	3.58	1.78	0.129	3.6	47.2	3.42	0.65	0.047	1.4

ACCURACY (CLSI EP9)⁽⁴⁾

The performance of this method (y) was compared with the performance of a similar salicylate method (x) on a Roche/Hitachi® 717. Thirty-nine patient serum samples ranging from 9.3-93.0 mg/dL (0.67-6.74 mmol/L) gave a correlation coefficient of 0.9986. Linear regression analysis gave the following equation:

$$\text{This method} = 1.02 (\text{reference method}) - 0.4 \text{ mg/dL} (0.03 \text{ mmol/L})$$

The performance of this method (y) with plasma was compared with the performance of this method (x) with serum on a Roche Hitachi® 717. Twenty-six patient serum and plasma samples spiked with salicylate ranging from 7.7-93.6 mg/dL (0.56-6.78 mmol/L) gave a correlation coefficient of 0.9972. Linear regression analysis gave the following equation:

$$\text{This method} = 1.01 [\text{Sekisui serum method}] - 0.7 \text{ mg/dL} (0.05 \text{ mmol/L})$$

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

The Americas
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

International
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

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Internationaal
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

Telefoon: 800-565-0265

Fax: 902-628-6504

E-mail: questions@sekisuidiagnostics.com

peidiagnostictechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

E-mail: info@sekisuidiagnostics.com

**Definitions for Symbols/ Définitions des symboles/
Definición de los símbolos/ Definizioni dei simboli/
Definitionen für Symbole/Beschrijving van symbolen**

This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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1. Done, A.K., *Salicylate Intoxication, Pediatrics*, 800-807 (1960).
2. Hammond, P.M., *A Simple Colorimetric Assay to Determine Salicylate Ingestion Utilizing Salicylate Mono-oxygenase*, VIIIth International Conference on Enzyme Engineering (1985).
3. You, K., Bittikofer, J.A., *Quantification of Salicylate in Serum by Use of Salicylate Hydroxylase*, Clin. Chem. 30/9, 1549-1551 (1984).
4. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
5. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, Third Edition (1990).
6. Kaplan, A., Szabo, L.L., *Clinical Chemistry: Interpretation and Techniques*, Lea and Febiger, Philadelphia (1983) p. 390.

Authorized Representative/Nom et adresse du représentant autorisé/Representante autorizado/Name e indirizzo del rappresentante autorizzato/Name and Adresse des zugelassenen Vertreters/Naam en adres van gevolmachtigde vertegenwoordiger:

Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF
United Kingdom
Tel (+44)(0)1732-220022
Fax (+44)(0)1732-220024

IN51120-10
June 7, 2011