## TOTAL BILIRUBIN L3K<sup>®</sup> ASSAY

**CATALOGUE NUMBER:** 295-10    **SIZE:** R1: 1 x 100 mL, R2: 1 x 25 mL  
295-30                                    R1: 3 x 100 mL, R2: 1 x 75 mL

### INTENDED USE

For the IN VITRO quantitative measurement of total bilirubin in serum and plasma.

### TEST SUMMARY

Bilirubin is a bile pigment which is normally found in serum as a result of red blood cell destruction. It is a product of hemoglobin breakdown by the reticuloendothelial system and exists in two forms. Unconjugated (indirect) bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and is excreted.

Elevated concentrations of total bilirubin may occur due to hemolytic processes, liver disease, or a disorder of the biliary tract.

Traditional methods of measuring serum bilirubin are based on the reaction of bilirubin with a diazo reagent to form the coloured compound azo-bilirubin. The diazo reaction can be accelerated by various chemicals. For example, Malloy-Evelyn<sup>(1)</sup> used ethanol, Jendrassik-Grof<sup>(2)</sup> used caffeine, and Walters-Gerarde<sup>(3)</sup> used DMSO. Modifications of these methods included the addition of surfactants as solubilising agents.<sup>(4)</sup>

### TEST PRINCIPLE

Bilirubin (both conjugated and unconjugated) couples with the diazo reagent in the presence of a surfactant to form azobilirubin. The increase of absorbance at 546 nm is directly proportional to the total bilirubin concentration.

### REAGENTS

Total Bilirubin L3K<sup>®</sup> Accelerator Reagent (R1): A solution containing buffer (pH of 1.1 at 25°C), 154 mmol/L NaCl, 190 mmol/L HCl, surfactants, and preservatives.

Total Bilirubin L3K<sup>®</sup> Diazo Reagent (R2): A solution containing buffer (pH of 0.9 at 25°C), 417 mmol/L HCl, 5 mmol/L 2,4 dichlorophenyldiazonium salt, and a surfactant.

### WARNINGS & PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.  
See Material Data Safety Sheet for additional information.

### REAGENT PREPARATION, STORAGE & STABILITY

The reagents are provided in a ready to use format.

The reagents are stable until the expiry date stated on the label when stored at 2-8°C. Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagents should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum, plasma and neo-natal serum samples. Sodium Heparin and Lithium Heparin specimens were tested and found acceptable. Use only suitable tubes or collection containers for specimen collection and storage. Various manufacturers of sample collection systems may contain differing materials that may affect test results. Follow all instructions provided by the collection system manufacturer.

### SAMPLE STORAGE

Specimens must be protected from light. Samples may be stored at 2-8°C for 3 days or at minus 70°C for 3 months in the dark.<sup>(5)</sup>

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(6)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from hemolysis, lipemia, and ascorbic acid were evaluated for this total bilirubin method on a Roche/Hitachi<sup>®</sup> analyzer using a significance criteria of >10% from control. Two concentrations of total bilirubin were evaluated. Interference data was collected in serum. Plasma data is expected to be similar.

Concentration of Total Bilirubin		Substance Tested	Concentration of Substance Where No Significant Interference was Observed	
mg/dL	µmol/L		mg/dL	µmol/L
1.4	23.9	Hemoglobin	1000	155
14.3	244.7		1000	155
1.4	24.1	Ascorbic Acid	µg/dL	µmol/L
13.3	226.9		3000	170
1.6	28.0	Intralipid	mg/dL	3000 mg/dL (33.9 mmol/L)
16.5	281.8		1000	Simulated Triglycerides

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(7)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

The reagents necessary for the determination of Bilirubin (Total) are provided.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material.
3. Quality control materials.

#### TEST CONDITION

For the data presented in this insert, studies using the Sekisui Diagnostics Total Bilirubin L3K<sup>®</sup> reagent were performed on an automated analyzer using a linear endpoint test mode with a sample to reagent ratio of 1:79 and a wavelength reading of 546 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration, if necessary, using an automated system is dependent on the system and parameters used.

#### QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

#### CALCULATIONS

The analyzer automatically calculates the total bilirubin concentration of each sample.

#### TEST LIMITATIONS

A sample with a total bilirubin concentration exceeding the linearity limit should be diluted with 0.9% saline and re-assayed incorporating the dilution factor in the calculation of the value.

#### REFERENCE INTERVALS<sup>(5)</sup>

0.2-1.0 mg/dL (3.4-17.1 µmol/L)

Neonates:	Premature	Full-Term
0-1 days	<8 mg/dL (<136.8 µmol/L)	2-6 mg/dL (34.2-102.6 µmol/L)
1-2 days	<12 mg/dL (<205.2 µmol/L)	6-10 mg/dL (102.6-171.0 µmol/L)
3-5 days	<14 mg/dL (<239.4 µmol/L)	4-8 mg/dL (68.4-136.8 µmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

#### PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi<sup>®</sup> analyzer unless otherwise stated.

#### RESULTS

Total Bilirubin concentration is reported as mg/dL (µmol/L).

#### REPORTABLE RANGE (CLSI EP6)<sup>(6)</sup>

The linearity of the procedure described is 35.0 mg/dL (598.5 µmol/L). The limit of quantitation of the procedure described is 0.1 mg/dL (1.7 µmol/L). This data results in a reportable range of 0.1 mg/dL to 35.0 mg/dL (1.7 – 598.5 µmol/L).

## PRECISION STUDIES (CLSI EP5)<sup>(6)</sup>

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Simple precision data was collected by assaying ten samples of two concentrations of control sera in 4 runs over 4 days.

Concentration		Total SD		Total	Concentration		Simple SD		Simple
mg/dL	µmol/L	mg/dL	µmol/L	CV%	mg/dL	µmol/L	mg/dL	µmol/L	CV%
0.9	15.4	0.05	0.9	5.4	0.9	15.4	0.02	0.3	2.5
4.6	78.7	0.21	3.6	4.6	4.5	77.0	0.11	1.9	2.5

Total precision data was collected on one concentration of control sera in 40 runs conducted over 10 days.

Within precision data was collected by assaying twenty samples of one concentration of control sera in one run.

Concentration		Total SD		Total	Concentration		Within Run SD		Within
mg/dL	µmol/L	mg/dL	µmol/L	CV%	mg/dL	µmol/L	mg/dL	µmol/L	Run
16.0	273.6	0.42	7.2	2.6	15.9	271.9	0.15	2.6	CV%
									1.0

## ACCURACY (CLSI EP9)<sup>(6)</sup>

The performance of this method (y) was compared with the performance of a similar method (x) on a Hitachi analyzer. Ninety-eight patient serum samples ranging from 0.2 to 32.0 mg/dL (3.4 to 547.2 µmol/L) gave a correlation coefficient of 0.9982. Linear regression analysis gave the following equation:

$$\text{This method} = 1.05 (\text{reference method}) + 0.3 \text{ mg/dL} (5.1 \text{ µmol/L}).$$

The performance of this method (plasma Na Heparin) was compared with the performance of a similar method (Serum) on a Hitachi analyzer. Thirty-nine patient serum and plasma samples ranging from 0.2 to 31.5 mg/dL (3.4 to 538.7 µmol/L) gave a correlation coefficient of 0.9851. Linear regression analysis gave the following equation:

$$\text{This method (plasma Na Heparin)} = 0.99 (\text{serum}) - 0.02 \text{ mg/dL} (0.3 \text{ µmol/L})$$

The performance of this method (plasma Lithium Heparin) was compared with the performance of a similar method (Serum) on a Hitachi analyzer. Forty-nine patient serum and plasma samples ranging from 0.1 to 34.3 mg/dL (1.7 to 586.5 µmol/L) gave a correlation coefficient of 0.9942. Linear regression analysis gave the following equation:

$$\text{This method (plasma Lithium Heparin)} = 1.06 (\text{serum}) - 0.2 \text{ mg/dL} (3.4 \text{ µmol/L})$$

The information stated above is based on results from Sekisui Diagnostics P.E.I. Inc. studies and is current at the date of publication.

## TRADEMARKS

L3K is a registered trademark of Sekisui. All other trademarks, brands, product names are the property of their respective companies

Manufactured by:

**The Americas**  
Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265  
Fax: 902-628-6504  
Email: questions@sekisuidiagnostics.com  
peidiagnostictechnical@sekisuidiagnostics.com

**International**  
Sekisui Diagnostics (UK) Limited  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, UK

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