UREA L3K[®] ASSAY

CATALOGUE NUMBER: 283-30 **SIZE:** 6 x 30 mL
 283-17 5 x 100 mL
 283-99 1 x 1000 mL

Note: Changes are highlighted.

INTENDED USE

For the IN VITRO quantitative measurement of urea in serum.



TEST SUMMARY

The measurement of urea can be clinically useful in the diagnosis of kidney dysfunction. Serum urea plays an important role in the discrimination between prerenal and postrenal azotemia.⁽¹⁾

The measurement of urea nitrogen has traditionally been performed by either condensation with diacetyl monoxime or by conversion of urea by urease to ammonia. Fearon⁽²⁾ first proposed the diacetyl monoxime method in 1939. This colorimetric method has limitations such as poor specificity and color instability. The use of urease in urea measurements was introduced by Marshall⁽³⁾ who measured the liberated ammonia by titration with acid. Ammonia produced by the urease action has also been measured by Nesslerization techniques^(4,5) and by the Berthelot reaction.⁽⁶⁾ These colorimetric methods lack specificity, require long incubation periods, and require high temperatures.

Talke and Schubert⁽⁷⁾ described the first totally enzymatic procedure. Urea L3K[®] is an enzymatic procedure employing glutamate dehydrogenase and a stabilized NADH analog,⁽⁸⁾ which is easy to use and applicable to routine instrumentation.

TEST PRINCIPLE

In the Talke and Schubert procedure, urea is first hydrolyzed by urease to give ammonia and carbon dioxide as in Equation I.

The ammonia produced in the first reaction reacts with 2-oxoglutarate and the reduced cofactor in the presence of glutamate dehydrogenase (GLDH) to yield glutamate and the cofactor (II). The decrease in the concentration of the reduced cofactor is monitored at 340 nm or 380 nm and is proportional to the concentration of the urea in the sample.

REAGENTS

Urea Reagent: A solution containing a buffer (pH 8.0 at 25°C), 14 mmol/L 2-oxoglutarate, 5.0 mmol/L ADP, > 12 KU/L GLDH (mammalian), > 50 KU/L Urease (botanical), 0.2 mmol/L NADH analog, stabilizers, and a preservative.

WARNINGS & PRECAUTIONS FOR USE

Avoid ingestion.
 S24/25: Avoid contact with skin and eyes.
 See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are provided in a ready to use format.

Supplied reagent is stable at 2-8°C until expiry date. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum. Do not use serum preserved with fluoride. Fluoride will inhibit the urease reaction.

SAMPLE STORAGE

Since urea is susceptible to bacterial degradation, specimens should be stored at 4-8°C until analysis.⁽⁹⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽¹⁰⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Urease is specific for urea, however, ammonia contamination will seriously affect the results obtained using this system. Analysis should not be performed in close proximity to a urinalysis laboratory or in a laboratory using cleaning supplies containing ammonia.

Interferences from icterus, lipemia, and hemolysis were evaluated for this urea method on a Roche/ Hitachi[®] 717 analyzer using a significance criterion of > 10% variance from control.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
mg/dL	mmol/L			
27	9.6	Hemoglobin	1000 mg/dL	155 µmol/L
29	10.3	Bilirubin	40 mg/dL	684 µmol/L
26	9.3	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽¹¹⁾

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' Urea L3K[®] reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material
3. Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a rate test mode, with a sample to reagent ratio of 1:100 and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the urea concentration of each sample.

TEST LIMITATIONS

A sample with a urea concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽⁹⁾

Urea
11-37 mg/dL (3.9-13.1 mmol/L)

Urea Nitrogen
5-17 mg/dL (1.8-6.0 mmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi[®] 717 analyzer unless otherwise stated.

RESULTS

Urea concentration is reported as mg/dL (mmol/L).

REPORTABLE RANGE (CLSI EP6)⁽¹⁰⁾

The linearity of the procedure described is 150 mg/dL (53.6 mmol/L). The lower limit of detection of the procedure described is 1 mg/dL (0.2 mmol/L). The limit of quantitation of the procedure described is 4 mg/dL (1.4 mmol/L). This data results in a reportable range of 4 to 150 mg/dL (1.4 to 53.6 mmol/L).

PRECISION STUDIES (CLSI EP5)⁽¹⁰⁾

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV%	Concentration		Within Run SD		Within Run CV%
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	mg/dL	mmol/L	
13.7	4.9	0.5	0.2	4.0	13.5	4.8	0.4	0.1	2.9
50.3	18.0	1.2	0.4	2.4	50.2	19.7	0.8	0.3	1.6

Within run precision data was collected on two concentrations of control sera each run 20 times in a single assay.

ACCURACY (CLSI EP9-P)⁽¹⁰⁾

The performance of this method (y) was compared with the performance of a similar urea method (x) on a Roche/Hitachi[®] 717 analyzer. Forty-two patient serum samples ranging from 4-61 mg/dL (1.4-21.8 mmol/L) were tested and gave a correlation coefficient of 0.9951. Linear regression analysis gave the following equation:

$$\text{This method} = 1.00 (\text{reference method}) - 0.1 \text{ mg/dL (0.04 mmol/L)}$$

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

SEKISUI
DIAGNOSTICS

The Americas
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
peidiagnosticttechnical@sekisuidiagnostics.com

International
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

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Internacional

Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, RU

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Temperature limitation
Límites de temperatura

REFERENCES/ REFERENCIAS

1. Tietz, N.W., *Textbook of Clinical Chemistry*, W.B. Saunders Co., 1268 (1986).
2. Fearon, W.R., *The Carbanido Diacetyl Reaction: A Test for Citrullin*, *Biochem. J.* 33:902 (1939).
3. Marshall, E.K., Jr., *A New Method for the Determination of Urea in Blood*, *J. Biol. Chem.* 15:487 (1913).
4. Gentzkow, C.J., *An Accurate Method for the Determination of Blood Urea Nitrogen by Direct Nesslerization*, *J. Biol. Chem.* 143:531 (1942).
5. Karr, W.B., *J. Lab. Clin. Med.* 9:329 (1924).
6. Fawcett, J.K., and Scott, J.E., *A Rapid and Precise Method for the Determination of Urea*, *J. Clin. Pathol.* 13:156 (1960).
7. Talke, H., Schubert, G.E., *Enzymatishche Harnstoffbestimmung in BLUT and Serum in Optischen Test NACH Warburg*, *Klin. Wchnschr* 43:174 (1965).
8. U.S. Patent No. 5,801,006.
9. Kaplan, L.A. and Pesce, A.J., *Clinical Chemistry - Theory, Analysis, and Correlation, Third Edition*. Mosby Year-Book Inc., St. Louis, p 500 (1996)
10. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
11. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, Third Edition (1990).

Authorized Representative/ Representante autorizado:
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF
United Kingdom
Tel (+44)(0)1732-220022
Fax (+44)(0)1732-220024

IN28330-10
July 15, 2011