Uric Acid-SL Assay

CATALOGUE NUMBER: 237-60 **SIZE:** 2 x 100 mL

INTENDED USE



For the IN VITRO quantitative measurement of uric acid in serum.

TEST SUMMARY

The measurement of serum uric acid is most helpful in the diagnosis of gout where the concentrations are chronically elevated. Hyperuricemia is also associated with lympho-proliferative disorders, chronic hemolytic anemias, and increased metabolism of nucleoproteins. Impaired renal function also leads to increased uric acid concentrations⁽¹⁾.

Uric acid is produced by the action of xanthine oxidase on xanthine and hypoxanthine which are products of nucleic acid degradation. The classic chemical method for the measurement of uric acid is based upon the reduction of phosphotungstic acid by uric acid to a blue phosphotungstate complex⁽²⁾. The method is not specific for uric acid because other reducing agents present in serum or plasma will also produce a blue phosphotungstate complex. In 1980 Fossati, et.al.⁽³⁾ described a procedure for assaying uric acid using uricase which produces hydrogen peroxide from uric acid. The hydrogen peroxide is then reacted with a phenolic compound to produce a red colored dye which can be measured spectrophotometrically in the visible range. The procedure is more specific and the phenolic compound, 3, 5-dichloro-2-hydroxybenzene sulfonate (DHBS), increases sensitivity due to the high absorbance coefficient of the quinoneimine dye produced. This uric acid procedure is a modification of that of Fossati⁽³⁾.

TEST PRINCIPLE

The hydrogen peroxide formed by the action of uricase on uric acid causes oxidative coupling of DHBS and 4-aminoantipyrine, in the presence of peroxidase, forming a red colored quinoneimine dye. The intensity of the color produced is directly proportional to the concentration of uric acid in the sample, with maximum absorbance at 520 nm.

REAGENTS

Uric Acid-SL Reagent: A buffered solution containing 1.8 mmol/L DHBS, 0.5 mmol/L 4-aminoantipyrine, > 3500 U/L peroxidase (botanical), > 200 U/L uricase (microbial), stabilizers, and preservatives.

WARNINGS & PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.
See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use.
Supplied reagent stable at 2 to 8°C until expiry date.
Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum.

SAMPLE STORAGE

Samples may be stored at 2 to 8°C for 3-5 days and for 6 months at minus 20-0°C⁽⁴⁾.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁶⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Lipemic samples will cause false elevation of uric acid values. Samples with visible lipemia should not be used with this procedure.

Interferences from icterus and hemolysis were evaluated for this Uric Acid method on a Roche/Hitachi[®] 704 analyzer using a significance criterion of >10% variance from control. Interference data was collected in serum.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
3.9 mg/dL	231 µmol/L	Hemoglobin	100 mg/dL	15.5 µmol/L
4.5 mg/dL	268 µmol/L	Bilirubin	8.0 mg/dL	136.8 µmol/L

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁵⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' Uric Acid-SL reagents.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

TEST CONDITION

For data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:75 and a wavelength reading (primary/secondary) of 505/660 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the uric acid concentration of each sample.

TEST LIMITATIONS

A sample with a uric acid concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽¹⁾

Males: 3.5-7.2 mg/dL (208-428 µmol/L)
Females: 2.6-6.0 mg/dL (155-357 µmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi[®] 704 automated analyzer unless otherwise stated.

RESULTS

Uric acid concentration is reported as mg/dL (µmol/L).

REPORTABLE RANGE (CLSI EP6)⁽⁶⁾

The linearity of the procedure described is 20.0 mg/dL (1190 µmol/L). The lower limit of detection is 0.3 mg/dL (18 µmol/L). This data results in a reportable range of 0.3-20.0 mg/dL (18-1190 µmol/L).

PRECISION STUDIES (CLSI EP5)⁽⁶⁾

Data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV %	Within Run SD		Within Run (or Simple) CV %
mg/dL	µmol/L	mg/dL	µmol/L		mg/dL	µmol/L	
4.13	246	0.1	6.0	2.5	0.05	3.0	1.3
7.43	442	0.2	11.9	2.0	0.07	4.2	1.0

ACCURACY (CLSI EP9)⁽⁶⁾

The performance of this method (y) was compared with the performance of a similar uric acid method (x) on a Roche/Hitachi[®] 704. Fifty patient serum samples ranging from 1.5 - 12.5 mg/dL (89 - 744 µmol/L) gave a correlation coefficient of 0.9937. Linear regression analysis gave the following equation:

$$\text{This method} = 1.078 (\text{reference method}) + 0.1 \text{ mg/dL (3 } \mu\text{mol/L)}.$$

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

The Americas
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

International
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

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