CREATINE KINASE-SL ASSAY

CATALOGUE NUMBER: 326-10 **SIZE:** R1: 1 x 100 mL, R2: 1 x 25 mL
326-30 R1: 3 x 100 mL, R2: 1 x 75 mL

INTENDED USE

For the IN VITRO quantitative measurement of creatine kinase activity in serum.

TEST SUMMARY

Elevated concentrations of serum creatine kinase are found in all types of muscular dystrophy and disease involving the destruction of muscle tissues. In the case of progressive muscular dystrophy, elevated CK concentrations may be present before clinical symptoms appear. Elevated CK concentrations are observed after myocardial infarction and in cases of acute cardiovascular disease and Reye's Syndrome. Neurogenic muscle diseases such as multiple sclerosis and polio do not produce elevated CK concentrations. CK activity is in inverse relation to thyroid activity, thus, in cases of hypothyroidism, elevated concentrations of CK would be observed.⁽¹⁾

The method of assaying CK using creatine phosphate and adenosine diphosphate (ADP) as substrates instead of creatine and adenosine triphosphate (ATP) was first described by Oliver.⁽²⁾ This method employs the optimized conditions as jointly developed by the Scandinavian Committee on Enzymes and the German Society for Clinical Chemistry.^(3,4,5) In the procedure N-acetylcysteine (NAC) is the thiol activator, and adenosine monophosphate (AMP) and p¹, p³-di(adenosine-5') pentaphosphate (Ap5A) are added to inhibit interference caused by adenylate kinase activity.

TEST PRINCIPLE

The reactions proceed as follows:

The rate of increase in absorbance at 340 nm due to the formation of NADPH is directly proportional to the creatine kinase activity.

REAGENTS

Creatine Kinase-SL Buffer Reagent (R1): A solution containing a buffer (pH 6.5), 27 mmol/L glucose, 27 mmol/L NAC, 14 mmol/L magnesium acetate, 2 mmol/L EDTA·Na₂, 2.7 mmol/L NADP, >5 KU/L hexokinase (yeast).

Creatine Kinase-SL Substrate Reagent (R2): A solution containing 11 mmol/L ADP, 28 mmol/L AMP, 55 μmol/L Ap5A, >14KU/L G-6-PDH (microbial), 2 mmol/L EDTA·Na₂, 160 mmol/L creatine phosphate.

WARNINGS & PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.
See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

The reagents are provided in a ready to use format. A working reagent may be prepared by mixing 4 parts CK-SL Buffer Reagent (R1) with 1 part CK-SL Substrate Reagent (R2).

The reagents included are stable until the expiry date stated on the labels at 2-8°C. Protect the reagent from sunlight. The working reagent is stable for 21 days at 2-8°C.

Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum.

SAMPLE STORAGE


The CK activity in serum decreases rapidly with increasing temperature and the sample should be separated and cooled to 2-8°C as quickly as possible. Macro, or mitochondrial CK, is the most stable isoenzyme followed by CK-3 (MM), CK-2 (MB) and CK-1 (BB). The average stabilities, defined as less than 5% loss of activity, are as follows:^(5, 1)

CK-MM (CK-3): 24 hours at 18-26°C, 1 week at 4°C, 1 month at -20°C
CK-MB (CK-2): 12 hours at 18-26°C, 3 days at 4°C, 1 month at -20°C
CK-BB (CK-1): less than CK-MB (CK-2)

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁶⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Adenylate kinase interferes with the determination of CK activity by generating additional ATP (see reaction below) and thus increasing the apparent creatine kinase activity.

The addition of AMP and Ap5A to the reagent in this procedure inhibits the adenylate kinase activity.⁽⁴⁾

Interference from icterus, lipemia and hemolysis were evaluated for this creatine kinase method on a Roche/Hitachi® analyzer using a significance criterion of > 10% variance from control. Interference data was collected in serum.

Concentration of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant	
198 U/L	Hemoglobin	200 mg/dL	31 μmol/L
205 U/L	Bilirubin	40 mg/dL	684 μmol/L
188 U/L	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

The information presented above is based on results from Genzyme Diagnostics studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁷⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Genzyme Diagnostics' Creatine Kinase-SL reagents.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
- 2) Calibration material. (If applicable.)
- 3) Quality Control materials.

TEST CONDITION

For data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:35 and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Genzyme Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

The frequency of calibration, if necessary, using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the creatine kinase activity of each sample.

TEST LIMITATIONS

A sample with a creatine kinase activity exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value. Note that dilution of highly active sera causes a progressive increase in activity per unit volume.⁽⁸⁾

REFERENCE INTERVALS⁽¹⁾

Males: 38-174 U/L (37°C)
Females: 26-140 U/L (37°C)

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® analyzer unless otherwise stated.

RESULTS

Creatine kinase activity is reported as U/L.

REPORTABLE RANGE (CLSI EP6)⁽⁶⁾

Reportable range using automated procedures will depend on the sample to reagent ratio used. The automated procedure described gives a reportable range from 2-1500 U/L.

PRECISION STUDIES (CLSI EP5)⁽⁶⁾

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Activity	Total SD	Total CV %	Within Run SD	Within Run CV %
U/L	U/L		U/L	
288	7.1	2.5	1.3	0.4
748	15.7	2.1	2.3	0.3

ACCURACY (CLSI EP9)⁽⁶⁾

The performance of this method (y) was compared with the performance of a similar method (x) on a Roche/Hitachi® analyzer. Fifty-five patient serum samples ranging from 31-1480 U/L gave a correlation coefficient of 0.9993. Linear regression analysis gave the following equation:

$$\text{This method} = 1.02 (\text{reference method}) + 1.3 \text{ U/L.}$$

The information presented above is based on results from Genzyme Diagnostics' studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

The Americas
Genzyme Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265
Fax: 902-628-6504
Email: customerservice@genzyme.com
peidiagnosticttechnical@genzyme.com

International
Genzyme Diagnostics
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

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