EN Genzyme Cystatin C Calibrator

INTENDED USE

For the calibration of the Genzyme Cystatin C assay.

SUMMARY AND PRINCIPLE

See the Genzyme Cystatin C Reagent package insert.

REAGENTS/COMPOSITION

The Genzyme Cystatin C calibrators consist of a bovine serum albumin liquid matrix and < 0.1% sodium azide with assigned concentrations of cystatin C.

The calibrators, labeled 1 through 6, have approximate values of 0.0, 0.5, 1.0, 2.0, 4.0 and 8.0 mg/L, respectively.

Precautions and Warnings

1. For in vitro diagnostic use only.
2. Do not use the calibrators after the expiration date printed on the labels.
3. **Caution:** The Genzyme Cystatin C calibrators contain < 0.1% sodium azide as an antimicrobial agent. Sodium azide may react with lead and copper plumbing to form potentially explosive metal azide buildup. Flush with copious amounts of water when discarding material.

Preparation of Calibrator

Calibrators: Liquid, ready to use.

Storage and Stability

Un-opened, the calibrators are stable until the expiration date shown on the label when stored at 2 - 8°C.

Once opened, and capped, the calibrators are stable for 12 weeks at 2 - 8°C.

DO NOT FREEZE.

Indications of Deterioration

Presence of turbidity or microbial growth may indicate deterioration.

PROCEDURE

Materials Provided

Description	Configuration	Catalog Number
Genzyme Cystatin C Calibrator	6 Levels x 1.5 mL vials	80-7149-00

Materials Required but not Provided

Description	Configuration	Catalog Number
Genzyme Cystatin C Reagent	R1 2 x 20 mL	80-7136-00
	R2 2 x 5 mL	
Genzyme Cystatin C Reagent	R1 2 x 48 mL	80-7143-00
	R2 2 x 12 mL	

- Quality Control materials.
- Analyzer capable of running two reagent chemistries.

Calibration

The values of the Genzyme Cystatin C calibrators can be found on the calibrator vial labels.

Refer to the instrument operator's manual for analyzer specific calibration procedures and for guidance in determining calibration frequency.

Quality control values should be within the expected ranges.

NOTE: To prevent potential contamination caused by intermixing the caps, uncap and dispense only one vial at a time and then recap immediately.

Manufactured by:

genzyme
Diagnostics

The Americas

Genzyme Corporation
31 New York Ave
Framingham, MA 01701-9322
USA
Phone: 800-332-1042
Fax: 800-762-6311

Email: usdiagcustomerservice@genzyme.com

International

Genzyme Diagnostics
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

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