LACTATE DEHYDROGENASE-SL ASSAY

CATALOGUE NUMBER: 327-10 **SIZE:** R1: 1 x 100 mL, R2: 1 x 25 mL
327-30 R1: 3 x 100 mL, R2: 1 x 75 mL

INTENDED USE


For the IN VITRO quantitative measurement of lactate dehydrogenase in serum.

TEST SUMMARY

Elevated concentrations of lactate dehydrogenase are clinically significant and can be found in disease states which result in cell damage. Myocardial infarcts, liver disease, megaloblastic anemias, renal disease, progressive muscular dystrophy and some malignancies all produce elevated lactate dehydrogenase values in serum.⁽¹⁾

Wacker et.al.⁽²⁾ published a method for the measurement of lactate dehydrogenase (LDH) utilizing lactate as the substrate and nicotinamide adenine dinucleotide (NAD) as the indicator coenzyme. The use of the reverse reaction, pyruvate to lactate (LDH-P), has been described by Wroblewski and LaDue.⁽³⁾ Amador et.al.⁽⁴⁾ claim the lactate to pyruvate (LDH-L) method is the method of choice because of the greater linearity of the reaction and improved stability of the reagents involved. The LDH-L reaction has been further studied by Gay et.al.⁽⁵⁾ and the optimum reaction conditions outlined. This procedure uses the LDH-L method of Wacker in accordance with the recommendations of the International Federation of Clinical Chemists.⁽⁶⁾

TEST PRINCIPLE

Lactate dehydrogenase converts L-lactate and NAD to pyruvate and NADH. The rate of increase in absorbance of the reaction mixture at 340 nm, due to the formation of NADH, is proportional to the LDH activity.

REAGENTS

LDH-SL Buffer Reagent (R1): A solution containing 325 mmol/L N-methyl-D-glucamine.

LDH-SL Substrate Reagent (R2): A solution containing 50 mmol/L L-lactate and 10 mmol/L nicotinamide adenine dinucleotide (NAD).

WARNINGS & PRECAUTIONS FOR USE

Avoid ingestion.

S24/25: Avoid contact with skin and eyes.

See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

The reagents are ready for use. A single working reagent can be prepared by mixing 4 parts of LDH-SL Buffer Reagent (R1) with 1 part of LDH-SL Substrate Reagent (R2).

Supplied reagent is stable until the expiry date when stored at 2-8°C and protected from light. The working reagent is stable for 24 hours at 2-8°C and for 4 hours at 18-26°C when protected from the light.

Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum. Specimens for analysis should be separated from the clot and assayed as soon as possible after collection.

SAMPLE STORAGE

Serum specimens stored at 18-26°C are stable for 2-3 days. Refrigerated storage may result in reduced LDH activity.⁽¹⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁷⁾

Cross Contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia and hemolysis were evaluated for this method on a Roche/Hitachi® 717 analyzer using significance criterion of > 10%. Interference data was collected in serum.

Red blood cells contain high concentrations of LDH, therefore, hemolyzed samples should not be used for LDH determination.

Activity of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant	
76 U/L	Bilirubin	40 mg/dL	684 µmol/L
71 U/L	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁸⁾

The information presented above is based on results from Genzyme Diagnostics' studies and is current at the date of publication.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Genzyme Diagnostics' LDH-SL reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration Material. (If applicable.)
3. Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:29 and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Genzyme Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

The frequency of calibration, if necessary, using an automated system is dependent on the system and parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the LDH activity of each sample.

TEST LIMITATIONS

A sample with a LDH activity exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution fraction in the calculation of the value.

REFERENCE INTERVALS⁽¹⁾

Adults: 100-210 U/L at 37°C

Higher concentrations are observed in children and infants.

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 704 analyzer unless otherwise stated.

RESULTS

LDH activity is reported as U/L.

REPORTABLE RANGE (CLSI EP6)⁽⁷⁾

The reportable range data as described was collected on a Roche/Hitachi® 717 analyzer. The linearity of the procedure described is 1350 U/L. The lower limit of detection of the procedure described is 10 U/L. This data results in a reportable range of 10 to 1350 U/L.

PRECISION STUDIES (CLSI EP5)⁽⁷⁾

Precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Activity (U/L)	Total SD (U/L)	Total CV %	Within Run SD (U/L)	Within Run CV %
173	2.3	1.3	1.6	0.9
433	6.0	1.4	2.6	0.6

ACCURACY (CLSI EP9)⁽⁷⁾

The performance of this method (y) was compared with the performance of a similar LDH method (x). Fifty-four patient serum samples ranging from 58- 2139 U/L were tested and gave a correlation coefficient of 0.9998. Linear regression analysis gave the following equation:

$$\text{This method} = 1.01 (\text{reference method}) - 4.8 \text{ U/L}$$

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

The Americas

Genzyme Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265

Fax: 902-628-6504

Email: customerservice@genzyme.com

peidiagnosticttechnical@genzyme.com

International

Genzyme Diagnostics
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

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Internacional

Genzyme Diagnostics
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, RU

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Temperature limitation
Límites de temperatura

REFERENCES/ REFERENCIAS

1. Burtis, C.A. and Ashwood, E.R. (Eds), Tietz Textbook of Clinical Chemistry, 2nd ed., W.B. Saunders Co., Philadelphia (1994).
2. Wacker, W.E.C., Ulmer, D.D., Vallee, B.L. Metalloenzymes and Myocardial Infarction, New England J. Med. 255, 449 (1956).
3. Wroblewski, F., LaDue, J.S., Lactate Dehydrogenase Activity in Blood, Proc. Soc. Exptl. Biol. Med. 90, 210 (1955).
4. Amador, E., Dorfman, L.E., Wacker W.E.C., Serum Lactic Dehydrogenase Activity: An Analytical Assessment of Current Assays, Clin. Chem. 9, 391-399 (1963).
5. Gay, R.J., McComb, R.B., Bowers, G.N., Optimum Reaction Conditions for Human Lactate Dehydrogenase Isoenzymes as They Affect Total Lactate Dehydrogenase Activity, Clin. Chem. 14, 740-753 (1968).
6. International Federation of Clinical Chemistry, Approved Recommendation on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes, Part 8. IFCC Method for Lactate Dehydrogenase. Eur. J. Clin. Chem. Clin. Biochem. 32, 639-655 (1994).
7. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.
8. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Third Edition, AACC Press, Washington, 1990.

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Genzyme Diagnostics
50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF
United Kingdom/ Reino Unido
Tel (+44)(0)1732-220022
Fax (+44)(0)1732-220024

IN32710-12
October 16, 2009