## MAGNESIUM (XB) ASSAY

**CATALOGUE NUMBER:** 125-12      **SIZE:** 4 x 125 mL  
125-50                                      1 x 1000 mL

**Note:** Changes are highlighted.

### INTENDED USE

For the IN VITRO quantitative measurement of magnesium in serum.

### TEST SUMMARY

Magnesium measurements are clinically important. Magnesium, along with potassium, is the most abundant intracellular cation. Magnesium is a required coenzyme in the metabolism of carbohydrate, lipids, and proteins. Decreased magnesium concentrations, hypomagnesemia, may cause tremors, muscle irritability, increased blood pressure and heart rate.

Elevated magnesium concentrations, hypermagnesemia, may lead to impaired respiration, coma, and cardiac arrest if left uncorrected.<sup>(1)</sup>

The method of choice for magnesium measurement is atomic absorption spectrophotometry.<sup>(2)</sup> This method is not practical for clinical laboratories so dye binding methods using methylthymol blue, titan yellow, calmagite, and xylidyl blue-1 are used instead. This method employs xylidyl blue-1 which allows the use of a single reagent which is rapid and precise.

### TEST PRINCIPLE

Xylidyl blue-1 + Mg<sup>++</sup> → Mg-xylidyl blue complex (red)

The absorbance increase, at 520 nm, of the red complex formed is directly proportional to the concentration of magnesium in the serum.

The absorbance decrease, at 660 nm, due to the complexing of the xylidyl blue-1 with Mg<sup>++</sup>, is directly proportional to the concentration of magnesium in the serum.

### REAGENTS

**Magnesium Reagent:** A solution containing buffer (pH 11.2 at 25°C), 0.14 mmol/L xylidyl blue-1, 0.1 mmol/L EGTA, and a surfactant.

**Magnesium Calibrator:** A solution containing 1.0 mmol/L (2.4 mg/dL) magnesium and a preservative (Not included with Cat. No. 125-50)

### WARNINGS & PRECAUTIONS FOR USE

S24/25: Avoid contact with skin and eyes.  
See Material Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use.

Supplied reagent is stable, unopened or opened, at 18-26°C until the expiry date stated on the label. Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum.

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(3)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

EGTA prevents calcium interference by preferential chelation of the calcium.

Interferences from icterus, lipemia, ascorbic acid, and hemolysis were evaluated for this magnesium method on a Roche/Hitachi® analyzer using a significance criterion of >10% variance from control. Interference data was collected in serum.

Concentration of Analyte		Substance Tested	Concentration of Interferent where Interference is Insignificant	
mg/dL	mmol/L			
1.9	0.78	Hemoglobin	800 mg/dL	124 µmol/L
1.9	0.78	Bilirubin	40 mg/dL	684 µmol/L
1.7	0.70	Ascorbic Acid	3000 µg/dL	170 µmol/L
2.0	0.82	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

The information presented above is based on results from Genzyme Diagnostics' studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(4)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Genzyme Diagnostics' Magnesium reagent.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

#### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:150, and wavelength readings of 520 or 600/660 nm. For assistance with applications on automated analyzers within Canada and the U.S, please contact Genzyme Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

#### QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

#### CALCULATIONS

The analyzer automatically calculates the magnesium concentration of each sample.

#### TEST LIMITATIONS

A sample with a magnesium concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed, incorporating the dilution factor in the calculation of the value.

#### REFERENCE INTERVALS<sup>(1)</sup>

1.6-2.6 mg/dL (0.65-1.05 mmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® analyzer unless otherwise stated.

### RESULTS

Magnesium concentration is reported as mg/dL (mmol/L).

To convert magnesium results to mEq/L, use the following conversion factor.<sup>(1)</sup>

$$\text{mmol/L} \times 2 = \text{mEq/L} \quad (\text{mEq/L} \times 0.5 = \text{mmol/L})$$

### REPORTABLE RANGE (CLSI EP6)<sup>(3)</sup>

The linearity of the procedure described is 6.1 mg/dL (2.5 mmol/L). The lower limit of detection of the procedure described is 0.5 mg/dL (0.2 mmol/L). This data results in a reportable range of 0.5 mg/dL (0.2 mmol/L) to 6.1 mg/dL (2.5 mmol/L). The reportable range will depend on the sample to reagent ratio used.

### PRECISION STUDIES (CLSI EP5)<sup>(3)</sup>

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV %	Concentration		Within Run SD		Within Run CV %
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	mg/dL	mmol/L	
2.0	0.82	0.1	0.05	5.9	1.8	0.75	0.1	0.03	4.5
4.0	1.65	0.2	0.06	3.8	4.0	1.64	0.1	0.04	2.4

Within run precision data was collected on two concentrations of control sera each run 20 times in a single assay.

### ACCURACY (CLSI EP9)<sup>(3)</sup>

The performance of this method (y) was compared with the performance of a similar xylidyl blue (x). Forty patient serum samples ranging from 0.9 mg/dL to 5.8 mg/dL (0.36-2.40 mmol/L) were tested and gave a correlation coefficient of 0.994. Linear regression analysis gave the following equation:

$$\text{This method} = 1.020 (\text{reference method}) - 0.2 \text{ mg/dL} (0.100 \text{ mmol/L}).$$

The information presented above is based on results from Genzyme Diagnostics studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**The Americas**  
Genzyme Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265  
Fax: 902-628-6504  
Email: customerservice@genzyme.com  
peidiagnosticttechnical@genzyme.com

**International**  
Genzyme Diagnostics  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, UK

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**Internacional**  
Genzyme Diagnostics  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, RU

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Manufacturer  
Fabricant  
Fabricante

Consult instructions for use  
Consulter les directives d'utilisation  
Consulte las instrucciones de uso

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Temperature limitation  
Limite de température  
Límites de temperatura

**REFERENCES/RÉFÉRENCES/ REFERENCIAS**

1. Tietz, N.W., (Editor) Clinical Guide to Laboratory Tests, W.B. Saunders Company, Philadelphia (1983) p.338.
2. Ratge, D., Kohse, K.P., and Wisser, H; Measurement of Magnesium in Serum and Urine with a Random Access Analyzer by use of a Modified Xylidyl Blue-1 Procedure, Clinical Chimica Acta, 159 (1986) 197-203.
3. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.
4. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, Third Edition, 1990.

IN12512-27  
June 3, 2010