



Maltase (α -D-Glucosidase)
Catalogue No. MALT-70-1235, 70-1235-01

Origin: Yeast


Specifications:

Specific Activity: > 50 U/mg protein @ 25°C
 Contaminants: (as α -Glucosidase activity = 100%)
 α - Galactosidase (< 0.01%)
 β - Galactosidase (< 0.01%)
 β - Glucosidase (< 0.01%)

Characteristics:


Molecular Weight:	~52kDa	
Optimum pH:	6-6.5 (K-PO ₄ buffer)	See Fig. 1
Optimum temp.:	35°C	See Fig. 3
pH stability:	6-8 (1 week, 5°C)	See Fig. 2
Thermal stability:	Below 30°C (pH7, 10 min)	See Fig. 4
Lyophilised stability:	2 years at -20°C	

Fig. 1 pH Optimum



● K-PO₄ ■ K-citrate □ TEA-HCL
 ○ Na- acetate △ Tris-HCL

Fig. 2 pH Stability




5°C, 1 week in 100mM buffer
 ● K-PO₄ ○ K-citrate ■ TEA-HCL
 □ GTA △ Tris-HCL

Please note that this product is for *in vitro* diagnostic use only


Maltase (α -Glucosidase)
(Catalogue No.1235)

Fig. 3 Temperature Optimum



75mM K-PO₄, pH 6.5, 138mM maltose, 5min

Fig. 4. Thermal Stability



10min, 0.1M K-PO₄, pH 7.0

K_m for Various Substrates:

Substrate	K _m Value
Maltose (G2)	$2.9 \times 10^{-2} \text{M}$
Maltotriose (G3)	$9.3 \times 10^{-3} \text{M}$
Maltotetraose (G4)	$9.5 \times 10^{-2} \text{M}$
Maltopentaose (G5)	$1.9 \times 10^{-1} \text{M}$


75 mM K-PO₄, pH 6.5, 25°C



Maltase (α -Glucosidase)
(Catalogue No. 1235)

Assay Principle:

Maltase catalyses the following reaction:



The generation of NADH can be measured spectrophotometrically at 340nm

Unit Definition:

One unit of activity is defined as the amount of enzyme which will catalyse the transformation of 1 μ mol of the substrate per minute under standard assay conditions.

(See Analytical Method for full details)

International
50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF, UK
Phone: +44 (0) 1732 220022
Fax: +44 (0) 1732 220024/5

The Americas
31 New York Avenue
Framingham, MA
01701-9322
Phone: 800 332 1042
Fax: 800 762 6311

Please note that this product is for *in vitro* diagnostic use only