EN

Sekisui Cystatin C Calibrator

INTENDED USE

For the calibration of the Sekisui Cystatin C assay.

SUMMARY AND PRINCIPLE

See the Sekisui Cystatin C Reagent package insert.

REAGENTS/COMPOSITION

The Sekisui Cystatin C calibrators consist of a bovine serum albumin liquid matrix and < 0.1% sodium azide with assigned concentrations of cystatin C.

The calibrators, labeled 1 through 6, have approximate values of 0.0, 0.5, 1.0, 2.0, 4.0 and 8.0 mg/L, respectively.

Precautions and Warnings

1. For *in vitro* diagnostic use only.
2. Do not use the calibrators after the expiration date printed on the labels.
3. **Caution:** The Sekisui Cystatin C calibrators contain < 0.1% sodium azide as an antimicrobial agent. Sodium azide may react with lead and copper plumbing to form potentially explosive metal azide buildup. Flush with copious amounts of water when discarding material.

Preparation of Calibrator

Calibrators: Liquid, ready to use.

Storage and Stability

Un-opened, the calibrators are stable until the expiration date shown on the label when stored at 2 - 8°C.

Once opened, and capped, the calibrators are stable for 12 weeks at 2 - 8°C.

DO NOT FREEZE.

Indications of Deterioration

Presence of turbidity or microbial growth may indicate deterioration.

PROCEDURE

Materials Provided

Description	Configuration	Catalog Number
Sekisui Cystatin C Calibrator	6 Levels x 1.5 mL vials	80-7149-00

Materials Required but not Provided

Description	Configuration	Catalog Number
Sekisui Cystatin C Reagent	R1 2 x 20 mL	80-7136-00
	R2 2 x 5 mL	
Sekisui Cystatin C Reagent	R1 2 x 48 mL	80-7143-00
	R2 2 x 12 mL	

- Quality Control materials.
- Analyzer capable of running two reagent chemistries.

Calibration

The values of the Sekisui Cystatin C calibrators can be found on the calibrator vial labels.

Refer to the instrument operator's manual for analyzer specific calibration procedures and for guidance in determining calibration frequency.

Quality control values should be within the expected ranges.


NOTE: To prevent potential contamination caused by intermixing the caps, uncap and dispense only one vial at a time and then recap immediately.

SEKISUI
DIAGNOSTICS

The Americas

Sekisui Diagnostics, LLC
31 New York Avenue
Framingham, MA 01701 USA
Phone: 800-332-1042
Fax: 800-762-6311
Email: questions@ sekisuidiagnostics.com

International

Sekisui Diagnostics (UK), Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK
Email: info@sekisuidiagnostics.com
www.sekisuidiagnostics.com

80-7163-00-01

5/12

©2012 Sekisui Diagnostics, LLC - All rights reserved.

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