OSOM® Influenza A&B Control Kit

Catolog number 191

INTENDED USE

The OSOM Influenza A&B Control Kit is intended for in vitro diagnostic use in quality control testing with the OSOM Influenza A&B Rapid Test.

SUMMARY

The OSOM Influenza A&B Control Kit provides swabs coated with a dried mixture containing influenza A or B antigens. The level of influenza A or B antigen applied to the swab is sufficient to provide a positive result when used with the OSOM Influenza A&B Rapid Test. The influenza A control swab acts as a negative control for the influenza B antigen and conversely, the influenza B control swab serves as a negative control for influenza A antigen. When used for quality control in the OSOM Influenza A&B Rapid Test, these controls are useful as an aid in verifying test and operator performance.

MATERIALS PROVIDED

5 Influenza A Positive Control Swabs (red shaft foam swab packaged with a desiccant tablet)

- Formalin inactivated influenza A/Kitakyushu/159/93 containing 0.05% sodium azide. Inactivity confirmed by inability of virus to infect cell culture.
- Result is representative of a mid-level positive.

5 Influenza B Positive Control Swabs (blue shaft foam swab packaged with a desiccant tablet)

- Formalin inactivated influenza B/Lee/40 containing 0.05% sodium azide. Inactivity confirmed by inability of virus to infect cell culture.
- Result is representative of a mid-level positive

MATERIALS REQUIRED BUT NOT PROVIDED

- OSOM Influenza A&B Rapid Test
- Clock or Timer

WARNINGS AND PRECAUTIONS

- DO NOT use the Controls beyond the expiration date.
- Potential biohazard. Handle as if potentially infectious. Handle this product according to established good laboratory practices using appropriate precautions.
- Not for patient use.

STORAGE

- The OSOM Influenza A&B Control Kit must be stored at room temperature (15°- 30°C; 59°- 86°F).

Stability

Swabs stored in the closed tube at 15°- 30°C (59°- 86°F) are stable until the expiration date printed on the tube and box label.

INSTRUCTIONS FOR USE

Use the OSOM Influenza A&B Control Swabs as you would a patient sample in accordance with the procedure for the OSOM Influenza A&B Rapid Test. Read the procedural instructions from the assay directional insert prior to testing.

EXPECTED RESULTS

The Positive A Control should yield a positive A result and a negative result for B. The Positive B control should yield a positive B result and a negative result for A. Refer to the assay directional insert for a full description in the Interpretation of Test Results section.

DISPOSING OF MATERIALS

Dispose of hazardous or biologically contaminated materials according to your institution’s practices. Discard all materials in a safe and acceptable manner that is in compliance with all country, state, and local requirements.

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