EN

## Albumin Assay

**CATALOGUE NUMBER:** 200-05      **SIZE:** 1 x 1000 mL  
 200-45                                      4 x 125 mL

### INTENDED USE

For the IN VITRO quantitative measurement of albumin in serum.

### TEST SUMMARY

Albumin's various functions in the blood stream include the transport of insoluble organic anions, the binding of toxic heavy metals, the transport of poorly soluble hormones, the maintenance of plasma colloidal osmotic pressure, and the provision of a reserve store of protein. Low serum values can result from malnutrition or liver disease, an increase in catabolism, increased excretion in urine or feces, or a change in distribution between the intravascular and extravascular compartments. High serum values can result from dehydration.<sup>(1)</sup>

Rodkey<sup>(2)</sup> in 1965 proposed a dye binding method for serum albumin in which he used 3',3'',5',5''-tetrabromo-m-cresolsulfonephthalein, sodium salt (Bromocresol green or BCG). Manual and automated methods for measuring serum albumin with BCG have been published by Hernandez, et.al.,<sup>(3)</sup> Dow and Pinto,<sup>(4)</sup> and Doumas et.al.<sup>(5)</sup> This procedure is a modification of that of Doumas exhibiting extended linearity.

### TEST PRINCIPLE

Albumin + BCG → BCG-Albumin complex

At an acid pH, the albumin binds to the BCG causing an increase in absorbance. The increase in absorbance at 630 nm due to the formation of the BCG-Albumin complex is directly proportional to the concentration of albumin in the sample.

### REAGENTS

Albumin Reagent: A solution containing acetate buffer (pH 4.2 at 25°C), ~ 0.39 mmol/L Bromocresol Green (BCG), a surfactant, and a preservative.

### WARNINGS & PRECAUTIONS FOR USE

Avoid contact with skin and eyes.

See Material Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use. Supplied reagent(s) stable at 18 to 26°C until expiry date. Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum.

### SAMPLE STORAGE

Specimens are stable at 18-26°C for up to one week, one month at 2-8°C, or for several months at 0 to -20°C.

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(6)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia, hemolysis, ascorbic acid and salicylate were evaluated for this albumin method on a Roche/Hitachi® 717 analyzer using a significance criterion of >10% from control.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
4.0 g/dL	40 g/L	Hemoglobin	600 mg/dL	93 µmol/L
3.8 g/dL	38 g/L	Bilirubin	40 mg/dL	684 µmol/L
4.2 g/dL	42 g/L	Ascorbic Acid	100 mg/dL	5678 µmol/L
3.8 g/dL	38 g/L	Salicylate	1 g/L	7.2 mmol/L
3.6 g/dL	36 g/L	Intralipid	800 mg/dL	2400 mg/dL (27.1 mmol/L) Simulated Triglycerides

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(7)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Sekisui Diagnostics' Albumin reagent.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material.
3. Quality Control materials.

#### TEST CONDITION

For data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:175 and a wavelength reading of (primary/secondary) 600/700 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

#### QUALITY CONTROL

A normal and abnormal concentration control serum should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

#### CALCULATIONS

The analyzer automatically calculates the albumin concentration of each sample.

#### TEST LIMITATIONS

A sample with an albumin concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

#### REFERENCE INTERVALS<sup>(8)</sup>

3.5-5.2 g/dL (35-52 g/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

#### PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 717 analyzer unless otherwise stated.

#### RESULTS

Albumin concentration is reported as g/dL (g/L).

## REPORTABLE RANGE (CLSI EP6)<sup>(6)</sup>

The linearity of the procedure described is 7.0 g/dL (70 g/L) in serum. The lower limit of detection is 0.0 g/dL (0 g/L). The limit of quantitation of the procedure described is 0.1 g/dL (1 g/L). This data results in a reportable range of 0.1 to 7.0 g/dL (1 to 70 g/L).

## PRECISION STUDIES (CLSI EP5)<sup>(6)</sup>

Data was collected on two control sera using a single lot of reagent in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV (%)	Concentration		Within SD		Within Run CV (%)
g/dL	g/L	g/dL	g/L		g/dL	g/L	g/dL	g/L	
3.9	39	0.08	0.8	2.1	3.9	39	0.06	0.6	1.4
2.7	27	0.05	0.5	1.7	2.7	27	0.04	0.4	1.5

Within run precision data was collected on two concentrations of control sera each run 20 times in a single assay.

## ACCURACY (CLSI EP9)<sup>(6)</sup>

The performance of this method (y) was compared with the performance of another albumin method (x) on a Roche/Hitachi® 717 analyzer. Forty patient serum samples ranging from 3.2 to 5.2 g/dL (32 to 52 g/L) were tested in duplicate and gave a correlation coefficient of 0.9965. Linear regression analysis gave the following equation:

$$\text{This method} = 1.04 (\text{reference method}) - 0.20 \text{ g/dL (2 g/L)}.$$

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**The Americas**  
Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265  
Fax: 902-628-6504  
Email: questions@sekisuidiagnostics.com  
peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

**International**  
Sekisui Diagnostics (UK) Limited  
Liphook Way  
Allington, Maidstone  
KENT, ME16 0LQ, UK

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The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.