EN

Creatinine-S Assay

CATALOGUE NUMBER: 221-30 SIZE: 2 x 250 mL + 1 x 125 mL
221-50 1 x 1000 mL + 1 x 250 mL

NOTE: Changes are highlighted.

INTENDED USE

For the IN VITRO quantitative measurement of creatinine in serum and urine.

TEST SUMMARY

In 1886 Jaffe developed an assay for creatinine based upon the reaction between creatinine and sodium picrate⁽¹⁾. In 1904 Folin⁽²⁾ used this reaction for the quantitative determination of creatinine in urine. Kinetic procedures based on the observed reaction rates of various substances, including creatinine, with alkaline picrate have been proposed by Fabing⁽³⁾ and Soldin⁽⁴⁾. This improved Jaffe chemistry is a kinetic procedure which does not require deproteinization of the sample and is formulated to reduce the interference of serum proteins.

Creatinine measurements are used in the diagnosis and treatment of renal diseases, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes.

TEST PRINCIPLE

Creatinine + alkaline picrate → creatinine-picric acid complex

At an alkaline pH, creatinine reacts with picric acid to form a Janovsky complex.

The rate of increase in absorbance at 510 nm due to the formation of the creatinine-picric acid complex is directly proportional to the concentration of creatinine in the sample.

REAGENTS

Creatinine Base Reagent (R1): A solution containing 0.25 mol/L sodium hydroxide and surfactants.

Creatinine Picric Acid Reagent (R2): A solution containing 31.4 mmol/L picric acid.

WARNINGS & PRECAUTIONS FOR USE

IVD

For *In Vitro* Diagnostic Use.

Rx ONLY

Warning

Contains: Sodium Hydroxide
H290 – May be corrosive to metals.
H319 – Causes serious eye irritation.
H315 – Causes skin irritation.

Prevention –

P264 – Wash thoroughly after handling.
P280 – Wear eye protection/face protection. **Wear protective gloves/protective clothing.**

Response – P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 – If eye irritation persists: Get medical advice/attention.
P302 + P350 – IF ON SKIN: Gently wash with plenty of soap and water.
P332 + P313 – If skin irritation occurs: Get medical advice/attention.
P362 – Take off contaminated clothing and wash before reuse.

Avoid contact with skin and eyes. The creatinine base reagent is caustic. The use of gloves while handling the reagent is recommended. See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use. R1 and R2 can be used separately on two reagent systems. A working reagent may be prepared by combining one volume of creatinine picric acid reagent and four volumes of creatinine base reagent. Mix well before using.

The reagents included are stable until the expiry date stated on the labels when stored at 18-26°C. The reagents may be stored in the refrigerated compartment of some analyzers. The working reagent is stable at 2-8°C for 18 months or the expiry stated on the labels, whichever is shorter. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration. For refrigerated conditions <5°C, precipitation may occur.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum or urine is the specimen of choice. Urine (24 hour) is to be collected without additives. For manual dilution, dilute urine with 0.9% saline or deionized water and multiply the result by the dilution factor. Pre-dilution of 1 in 10 to 1 in 20 is usually required for urine samples.

SAMPLE STORAGE

Specimens may be stored at 2-8°C for three days and at minus 20-0°C for longer periods.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁵⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia, hemolysis were evaluated for this creatinine method on a Roche/Hitachi® 704 analyzer using a significance criterion of > 10% variance from control. Interference data was collected in serum.

| Concentration of Analyte | | Substance Tested | Concentration of Interferent Where Interference is Insignificant | |
|--------------------------|-----------|------------------|--|--|
| Conventional Units | SI Units | | | |
| 1.0 mg/dL | 88 µmol/L | Hemoglobin | 750 mg/dL | 116 µmol/L |
| 1.0 mg/dL | 88 µmol/L | Bilirubin | 10 mg/dL | 171 µmol/L |
| 1.0 mg/dL | 88 µmol/L | Intralipid | 1000 mg/dL | 3000 mg/dL (33.9 mmol/L) Simulated Triglycerides |

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁶⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' creatinine reagents.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:24, and a wavelength reading of 505 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control serum or urine control should be analyzed as required. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the creatinine concentration of each sample.

TEST LIMITATIONS

A sample with a creatinine concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽⁷⁾

0.5-1.2 mg/dL (44-106 µmol/L)

Urine 24 Hour: Male: 800-2000 mg/24 hours (7072-17680 µmol/24 hours)
Female: 600-1800 mg/24 hours (5304-15912 µmol/24 hours)

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

PERFORMANCE CHARACTERISTICS

Data presented was collected on an automated analyzer unless otherwise stated.

RESULTS

Creatinine concentration is reported as mg/dL (µmol/L).

REPORTABLE RANGE (CLSI EP6)⁽⁵⁾

The linearity of the procedure described is 22.0 mg/dL (1945 µmol/L). The lower limit of detection of the procedure described is 0.1 mg/dL (9 µmol/L). This data results in a reportable range of 0.1-22.0 mg/dL (9-1945 µmol/L).

PRECISION STUDIES (CLSI EP5)⁽⁵⁾

Day to day precision was established by assaying two serum controls and two urine controls twice a day for 20 days.

| Creatinine-S | N | Mean | | Standard Deviation | | Coefficient of Variation % |
|--------------|----|-------|--------|--------------------|--------|----------------------------|
| | | mg/dL | µmol/L | mg/dL | µmol/L | |
| Serum 1 | 40 | 1.0 | 88 | 0.03 | 2.6 | 2.9 |
| Serum 2 | 40 | 6.6 | 583 | 0.15 | 13.3 | 2.3 |
| Urine 1 | 40 | 78.3 | 6922 | 1.09 | 96.4 | 1.4 |
| Urine 2 | 40 | 212.6 | 18,793 | 1.88 | 166.2 | 0.9 |

Within run precision was established by assaying two serum controls and two urine controls 20 times.

| Creatinine-S | N | Mean | | Standard Deviation | | Coefficient of Variation % |
|--------------|----|-------|--------|--------------------|--------|----------------------------|
| | | mg/dL | µmol/L | mg/dL | µmol/L | |
| Serum 1 | 20 | 1.0 | 88 | 0.02 | 1.8 | 1.7 |
| Serum 2 | 20 | 6.5 | 575 | 0.04 | 3.5 | 0.6 |
| Urine 1 | 20 | 78.3 | 6922 | 0.73 | 64.5 | 0.9 |
| Urine 2 | 20 | 212.6 | 18,793 | 1.17 | 103.4 | 0.5 |

ACCURACY (CLSI EP9)⁽⁵⁾

A comparison was made between this method and a similar creatinine procedure using 40 serum samples ranging from 0.74-20.17 mg/dL (65-1783 µmol/L). The correlation coefficient was 0.999. Linear regression analysis gave the following equation:

$$\text{This method} = 1.008 (\text{comparison method}) + 0.003 \text{ mg/dL (0.27 } \mu\text{mol/L)}.$$

A comparison was made between this method and a similar creatinine procedure using 40 urine samples ranging from 0.60-15.10 mg/dL (53-1335 µmol/L). The correlation coefficient was 0.9998. Linear regression analysis gave the following equation:

$$\text{This method} = 0.9535 (\text{comparison method}) + 0.01 \text{ mg/dL (0.88 } \mu\text{mol/L)}.$$

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

SEKISUI
DIAGNOSTICS

The Americas
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

International
Sekisui Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK

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The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.
MDSS GmbH