

ENZYMES

Reverse Transcriptase

ORIGIN *E.coli* (Recombinant)

CAT# TRT Range
EC# 2.7.7.7

SPECIFICATION

Concentration/Activity* 100 U/ μ L
RNase Activity None detected

*One unit is defined as the amount of enzyme required to incorporate 1 nmole of dTTP into an acid-insoluble material in 10 min at 42°C.

PRODUCT FORMAT

Formulated in 50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 0.1 mM EDTA, 10 mM DTT, 0.01% Nonidet P-40, 50% Glycerol and store at -20°C.

PART #	DESCRIPTION	CONTENT/FORMAT
TRT-101	ReverTra Ace RTase Kit	0.1mL (10,000U) Reverse Transcriptase 1mL 5x RTase Buffer
<i>On Request</i>	ReverTra Ace RTase 50KU	0.5mL (50,000U) Reverse Transcriptase
TRT-1B	5x RTase Buffer 1mL	1mL 5x RTase Buffer

DESCRIPTION AND APPLICATION


ReverTra Ace is a high efficient M-MLV (Moloney Murine Leukemia Virus) reverse transcriptase that has been genetically modified to remove RNase H activity and increase reaction efficiency. It is the preferred enzyme for applications requiring full-length cDNAs and high product yields from total RNA, mRNA, rRNA, etc.

CHARACTERISTICS


Features:

- RNase minus M-MLV RTase with improved performance.
- Enables the synthesis of longer cDNAs (\geq 14 kb) than the WT-enzyme.
- Exhibits excellent reaction efficiency at high temperatures.

Application Data: cDNAs were synthesized with oligo (dT)30 primers and 100 U enzyme/poly (A)+ RNA mixture (1.35-9.49 kb, 0.4 mg) as templates for 30 min at various temperatures. cDNAs were labeled with (32p-dCTP) during the reaction. The synthesized cDNAs were separated by 1% denatured agarose gel electrophoresis, and detected. The results suggested that our Reverse Transcriptase can elongate efficiently at 42-55°C compared to other RNase H minus RTases from other companies.



Application Data: G3PDH genes (500 bp) were amplified by PCR using cDNA templates that were synthesized with various RNase H minus RTases from G3PDH mRNA (10²-10⁵ copies/reaction). The RTase reaction was performed with specific reverse primers and 100 U enzyme at 42°C for 20 min. The results indicated that our Reverse Transcriptase is suitable for RT-PCR amplifications that require sensitivity.



Application Data: cDNA was synthesized by ReverTra Ace using a specific primer for the 3'-end of dystrophin mRNA at 42°C for 30 min. The 5' region at a distance of 14 kb from the 3' end of the dystrophin gene was amplified by PCR. The result indicated that ReverTra Ace can elongate cDNA of ≥ 14 kb.



M: Fx174/Hinc II Marker
1: ReverTra Ace
2: Company A
3: Company B

THE AMERICAS
 Sekisui Diagnostics, LLC
 4 Hartwell Place
 Lexington, MA 02421
 Phone: 800 332 1042
 Fax: 800 762 6311

INTERNATIONAL
 Sekisui Diagnostics (UK) Limited
 Liphook Way, Allington
 Maidstone, Kent, ME16 0LQ, UK
 Phone: +44 1622 607800
 Fax: +44 1622 607801

MANUFACTURED BY

DISTRIBUTED BY

Because every result matters™

engage@sekisuienzymes.com
 www.sekisuidiagnostics.com