## ACETAMINOPHEN L3K® ASSAY

**CATALOGUE NUMBER:** 506-10 **SIZE:** R1: 1 x 10 mL, R2: 2 x 10 mL  
506-30 R1: 3 x 10 mL, R2: 6 x 10 mL

**Note:** Changes are highlighted.

### INTENDED USE

For the IN VITRO quantitative measurement of acetaminophen in serum and plasma. Measurement of acetaminophen is used in the diagnosis and treatment of acetaminophen overdose toxicity.

### TEST SUMMARY


Acetaminophen (paracetamol) is used as an analgesic in many different formulations<sup>(1)</sup>. While therapeutic doses rarely cause adverse side effects, the effect of long term treatment with acetaminophen is unclear. Cases have been reported where chronic excessive use of acetaminophen has led to hepatotoxicity and nephrotoxicity.<sup>(2,3)</sup> In cases of acute overdosage, acetaminophen can cause severe hepatic damage leading to hepatic failure if untreated.<sup>(4,5,6)</sup>

The management of acetaminophen overdose requires early recognition of the drug in the bloodstream. Toxicity is generally reported at concentrations over 200 µg/mL (1324 µmol/L). N-acetylcysteine has been used as an antidote in conjunction with intensive support care. Early diagnosis of acetaminophen-induced hepatotoxicity is important since initiation of therapy within 8 hours of ingestion lessens the potential for hepatic injury, and decreases the mortality rate.<sup>(7)</sup>

The majority of methods for measuring acetaminophen are based on spectrophotometric or chromatographic principles. Chromatographic methods are specific for the parent compound, however, they are not well suited to emergency laboratories. Spectrophotometric methods are simpler and more rapid, but do not always offer the desired specificity.

This spectrophotometric method is rapid, reliable, convenient, and specific for acetaminophen.

### TEST PRINCIPLE

The enzyme, acyl amidohydrolase, cleaves the amide bond of the acetaminophen molecule, leaving p-aminophenol and acetate. The p-aminophenol is reacted with 2,5-dimethylphenol in the presence of manganese ions to form a colored compound, 4-(4-aminophenyl)-2,5-dimethylcyclohexadiene-1-one. The increased absorbance at 605 nm due to the formation of 4-(4-aminophenyl)-2,5-dimethylcyclohexadiene-1-one is directly proportional to the concentration of acetaminophen in the sample.

### REAGENTS

**Acetaminophen Enzyme Reagent (R1):** A solution containing buffer (pH 8.6 at 25°C), 0.3 mmol/L MnCl<sub>2</sub>·4H<sub>2</sub>O, ≥ 0.9 KU/L Acyl Amidohydrolase (microbial), 50 mg/L sodium azide.

**Acetaminophen Color Reagent (R2):** A solution containing 0.1 mol/L sodium carbonate buffer (pH 11.5 at 25°C), 61 mmol/L 2,5-dimethylphenol, stabilizer, preservative.

**Acetaminophen Calibrator:** 1 x 5 mL of a solution containing buffer (pH 5.2 at 25°C), 151 µg/mL (1000 µmol/L) acetaminophen, preservatives.

Internal reference standards are created for Acetaminophen using a USP grade reference Acetaminophen material (not less than 98% and not more than 102% of paracetamol on an anhydrous basis). Acetaminophen calibrator is manufactured gravimetrically and tested against these internal reference standards.

### WARNINGS AND PRECAUTIONS FOR USE

#### DANGER

Contains: Sodium Hydroxide

H314 – Causes severe skin burns and eye damage.

Prevention – P260 – Do not breathe vapor/spray.

P264 – Wash thoroughly after handling.

P280 – Wear protective gloves/protective clothing/eye protection/ face protection.

Response – P301 + P330 + P331 – IF SWALLOWED; rinse mouth. Do NOT induce vomiting.

P303 + P361 + P353 – IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P363 – Wash contaminated clothing before reuse.

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 – Immediately call a POISON CENTRE or doctor/physician.

Storage – P405 – Store locked up.

Disposal – P501 – Dispose of contents/container in accordance with local/regional /national/international regulations.

See Material Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use.

Supplied reagents are stable at 2-8°C until expiry date. Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagents should be clear. Turbidity would indicate deterioration. The presence of crystals would indicate deterioration; product with crystals should be replaced with fresh product.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum or lithium heparinized plasma. Separated samples may be stored for up to 14 days at 4 to 8°C prior to being tested. If testing will be delayed more than 14 days, separated samples may be stored frozen at ≤ -20°C for up to 45 days.<sup>(8)</sup>

### LIMITATIONS/ INTERFERING SUBSTANCES (CLSI EP7)<sup>(9)</sup>

Interferences from hemolysis, icterus and lipemia were evaluated for this acetaminophen method on a Roche/Hitachi® 717 using a significance criterion of > 10% variance from control.

Hemoglobin concentration greater than 200 mg/dL (31 µmol/L) showed a positive bias of up to 45% at acetaminophen concentration of 14.0 µg/mL (93 µmol/L). Hemoglobin produces significant interference in this method; therefore hemolyzed samples should not be used.

Intralipid concentration greater than 200 mg/dL showed a positive bias of up to 38% at acetaminophen concentration of 15.3 µg/mL (101 µmol/L). Lipemic samples should not be used.

Conjugated bilirubin concentration of up to 2 mg/dL (34.2 µmol/L) did not interfere in samples with acetaminophen concentrations of 16.3 µg/mL (108 µmol/L). Unconjugated bilirubin concentration of up to 2 mg/dL (34.2 µmol/L) did not interfere in samples with acetaminophen concentrations of 16.3 µg/mL (108 µmol/L).

**NOTE:** Significantly reduced acetaminophen recovery has been demonstrated in situations where testing for acetaminophen toxicity has been performed on hyperbilirubinemic samples at acetaminophen levels in the range of 15.1 mg/L (100 µmol/L). This interference is not detected at acetaminophen levels in the range of 45.3 mg/L (300 µmol/L) or higher. It is recommended that laboratories review the Rumack-Matthews Nomogram for patient ingestion status, treatment and monitoring protocols to determine the extent of the interference.

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(9)</sup>

Cross Contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia, hemolysis, ascorbic acid and N-acetylcysteine were evaluated for this acetaminophen method on the Roche/Hitachi® 717 analyzer using a significance criterion of >10% or ±1.25 µg/mL (8 µmol/L) variance from control, whichever is greater. Plasma data is expected to be similar.

Substance Tested	Concentration with no Significant Interference	Acetaminophen level
Hemoglobin	200 mg/dL (31 µmol/L)*	14.0 µg/mL (93 µmol/L)*
Conjugated Bilirubin	2 mg/dL (34.2 µmol/L)*	16.3 µg/mL (108 µmol/L)*
Unconjugated Bilirubin	2 mg/dL (34.2 µmol/L)	16.3 µg/mL (108 µmol/L)
Ascorbic Acid	3000 µg/dL (170 µmol/L)	15.7 µg/mL (104 µmol/L)
N-Acetylcysteine	1500 mg/L (9.2 mmol/L)	14.7 µg/mL (97 µmol/L)
Intralipid	200 mg/dL [600 mg/dL (6.8 mmol/L) Simulated Triglycerides]*	15.3 µg/mL (101 µmol/L)*

\* See additional information under the heading "Limitations/ Interfering Substances".

### ANALYTICAL SPECIFICITY TO DRUGS

Interferences from the following therapeutic drugs were tested at acetaminophen concentrations of 5.0 µg/mL (33 µmol/L) and 30.0 µg/mL (199 µmol/L) and were evaluated for this Acetaminophen method on a Roche/Hitachi® 717 analyzer using a significance criterion of >10% or ±1.25 µg/mL (8 µmol/L) variance from control, whichever is greater.

Substance Tested	Concentration with No Significant Interference
Theophylline	222 µmol/L
Phenylbutazone	2.89 mmol/L
Ibuprofen	2425 µmol/L
Imipramine	2.5 µmol/L
Acetylsalicylic Acid	6.51 mmol/L
Levodopa	25.3 µmol/L
Ampicillin	152 µmol/L
Doxycycline	67.5 µmol/L
Amitriptyline	3.61 µmol/L
Metronidazole	701 µmol/L
Cefoxitin	1546 µmol/L
Cyclosporin	10.0 µmol/L
Methyl-l-Dopa	71 µmol/L
Rifampicin	78.1 µmol/L
Salicylate	4.34 mmol/L
Ascorbic Acid	342 µmol/L

Samples containing NAPQ1 (N-Acetyl-4-benzoquinoneimine) may cause elevated levels of measured acetaminophen. Samples containing ≥20 mg/L metazolone may cause elevated levels of measured acetaminophen.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(10)</sup>

### ANALYTICAL PROCEDURE

#### MATERIAL PROVIDED

Sekisui Diagnostics' Acetaminophen reagents and calibrator.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

- Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
- Acetaminophen L3K Reagent Application files for Automated Analyzer.
- Quality control materials.

#### TEST CONDITIONS

For data presented in this insert, studies using Sekisui Diagnostics acetaminophen reagent were performed on an

automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:41 and a wavelength reading of 660 nm.

Beckman AU680 analyzer has been validated for the use with the Acetaminophen L3K Assay. Please refer to the application information.

For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

An acetaminophen calibrator is included and should be used as directed to calibrate the procedure. The frequency of calibration on automated systems is dependent on the system and the parameters used.

#### QUALITY CONTROL

Appropriate concentrations of quality control materials should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

#### CALCULATIONS

The analyzer calculates the acetaminophen concentration of each sample.

#### TEST LIMITATIONS

A sample with an acetaminophen concentration exceeding the linearity limit should be diluted with 0.9% saline and re-assayed incorporating the dilution factor into the calculation of the value.

#### REFERENCE INTERVALS<sup>(7)</sup>

Therapeutic concentration: 10-30 µg/mL (66-199 µmol/L)  
Toxic concentration: > 200 µg/mL (1324 µmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

#### PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 717 analyzer unless otherwise stated.

#### RESULTS

Acetaminophen concentration is reported as µg/mL (µmol/L).

To convert acetaminophen results to mg/L (µg/mL) or mg/dL, use the following conversion factors:  
µmol/L x 0.151 = mg/L (µg/mL)  
mg/dL x 10 = mg/L (µg/mL)

NOTE: 1 mg/L = 1 µg/mL

#### REPORTABLE RANGE (CLSI EP6)<sup>(9)</sup>

The linearity of the procedure described is 377.5 µg/mL (2500 µmol/L). The limit of quantitation of the procedure described is 0.6 µg/mL (4 µmol/L). This data results in a reportable range of 0.6 to 377.5 µg/mL (4 to 2500 µmol/L).

#### PRECISION STUDIES (CLSI EP5)<sup>(9)</sup>

Total precision data was collected on three control sera using a single lot of reagent in 40 runs conducted over 20 days. Within run precision data was collected by assaying twenty samples of three concentrations of control sera in one run using one lot of reagent.

Concentration		Total SD		Total CV %	Concentration		Within Run SD		Within Run CV%
µg/mL	µmol/L	µg/mL	µmol/L		µg/mL	µmol/L	µg/mL	µmol/L	
10.1	67	0.29	1.9	2.9	10.1	67	0.15	1.0	1.5
36.7	243	0.47	3.1	1.3	36.2	240	0.29	1.9	0.8
112.2	743	1.49	9.9	1.3	110.1	729	0.69	4.6	0.6

Additional precision analysis was conducted on two elevated concentrations of acetaminophen in sera. Total precision was collected over a 10 day period with 4 runs per day with each concentration done in duplicate.

Concentration		Total SD		Total CV%
µg/mL	µmol/L	µg/mL	µmol/L	
201.0	1331	2.6	17	1.3
320.2	2120	4.7	31	1.4

#### ACCURACY (CLSI EP9)<sup>(9)</sup>

The performance of this method (y) was compared with the performance of a similar acetaminophen method (x) on a Roche/Hitachi® 717. A combination of eighty-eight natural and spiked patient serum samples ranging from 5.7-356.5 µg/mL (38-2361 µmol/L) gave a correlation coefficient of 0.9998. Linear regression analysis gave the following equation:

$$\text{This method} = 1.064 (\text{reference method}) + 1.1 \mu\text{g/mL} (7.0 \mu\text{mol/L}).$$

The performance of this method with plasma (y) was compared to the performance of this method with serum (x) on a Roche/Hitachi® 717. Twenty-five serum and plasma samples spiked with acetaminophen ranging from 4.5 to 368.6 µg/mL (30 to 2441 µmol/L) gave a correlation coefficient of 0.9999. Linear regression analysis gave the following equation:

$$\text{This method (plasma)} = 0.999 [\text{This method (serum)}] - 0.3 \mu\text{g/mL} (2.2 \mu\text{mol/L})$$

#### TRADEMARK

L3K is a registered trademark of Sekisui. All other trademarks, brands, product names are the property of their respective companies.

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The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.
Sekisui Diagnostics (UK) Limited