# EN

## Sekisui Cystatin C Reagent

### INTENDED USE

For the quantitative measurement of cystatin C concentration in human serum, heparinized plasma, or EDTA plasma. Cystatin C measurements are used as an aid in the diagnosis and treatment of renal disease.

### SUMMARY

Cystatin C is a non-glycosylated 13 kDa basic protein that is a member of the cystatin super-family of cysteine proteinase inhibitors<sup>1-4</sup>. It is produced by most nucleated cells in a consistent manner, and unlike creatinine its rate of production is not affected by factors such as muscle mass or nutrition<sup>5</sup>. Due to its low molecular weight it is freely filtered by the healthy glomerulus and is neither secreted nor taken up by cells lining the nephron. This, together with its stable production rate, means that its serum concentration is strongly related to the glomerular filtration rate (GFR). In addition, its response to a change in GFR is sharper than that of creatinine<sup>6-7</sup>. Therefore, cystatin C is considered a useful endogenous marker for the detection of mild renal impairment<sup>6-7</sup>.

### PRINCIPLE

The Sekisui Cystatin C reagent is based on the sol particle immunoassay principle<sup>8-10</sup>. It contains colloidal gold particles coated with anti-cystatin C specific polyclonal antibodies. The reaction between the particles and any cystatin C in the sample results in the formation of agglutinates and an associated change in absorbance signal. The change in absorbance signal is proportional to the amount of cystatin C in the sample. Cystatin C concentration in the sample is determined by comparison with a standard curve.

### REAGENTS

#### Composition

Component	Ingredients	Concentration
Reagent 1	Buffer Solution Sodium Azide	pH 6.7 <0.1%
Reagent 2	Colloidal gold coated with rabbit polyclonal antibodies. Sodium Azide	<0.1%

#### Precautions and Warnings

1. For *in vitro* diagnostic use only.
2. Do not use the reagents beyond the expiration date printed on the label.
3. **Warning:** All specimens used in the test should be considered potentially infectious. Universal precautions, as they apply to your facility, should be adhered to for the handling and disposal of materials during and after testing<sup>11</sup>.
4. Reagents must be used with the Sekisui Cystatin C Calibrator.
5. **Caution:** Avoid freezing reagents
6. **Caution:** Reagent 1 and 2 contain <0.1% sodium azide as an antimicrobial agent. Sodium azide may react with lead and copper plumbing to form a potentially explosive

metal azide build-up. Flush with copious amounts of water when discarding material.

#### Reagent Preparation

Reagent 1: Liquid. Mix contents gently, prior to use.

Reagent 2: Liquid. Mix thoroughly by inverting the bottle a minimum of 20 times until the colloidal gold particles become a homogeneous suspension. Invert daily.

#### Storage and Stability

Unopened reagent is stable until the expiration date shown on the label when stored at 2 - 8°C.

Once opened, the reagent is stable for up to one month at 2 - 8°C.

DO NOT FREEZE.

#### Indications of Deterioration

Presence of turbidity or microbial growth may indicate deterioration.

Inability to recover control values.

### SPECIMEN COLLECTION AND PREPARATION

Serum, EDTA plasma, or lithium heparinized plasma are the recommended sample types. Use standard sample collection and preparation methods<sup>12</sup>.


If not analyzed promptly, serum or plasma specimens may be stored at 2-8°C for 19 days, or at 20 - 25°C for 4 days. If storage periods will be longer, the specimen should be frozen (-20°C).

Samples may be frozen and thawed twice.

### PROCEDURE

#### Assay

Below is a general example of the Sekisui Cystatin C assay procedure for an automated analyzer. All analyzer applications should be validated.

For assistance with applications on automated analyzers inside the US, please contact Sekisui Diagnostics Technical Marketing at (800) 332-1042. Outside the US, please contact your local distributor.

#### Materials Provided

Sekisui Diagnostics Cystatin C reagents 1 and 2 are required for the measurement of cystatin C.

Description	Configuration	Catalog Number
Sekisui Cystatin C Reagent	R1 2 x 20 mL	80-7136-00
	R2 2 x 5 mL	
Sekisui Cystatin C Reagent	R1 2 x 48 mL	80-7143-00
	R2 2 x 12 mL	

## Materials Required but not Provided

Description	Configuration	Catalog Number
Sekisui Cystatin C Calibrator	6 Levels x 1.5 mL vials	80-7149-00

- Quality Control materials.
- Analyzer capable of running two reagent chemistry.

### Calibration

Only the Sekisui Cystatin C Calibrator should be used to calibrate the Sekisui Cystatin C assay. The calibration curve is stable for at least 14 days with daily inversion of Reagent 2.

Refer to the instrument operator's manual for analyzer specific calibration procedures and for guidance on determining calibration frequency.

### Quality Control

Reliability of test results should be monitored routinely with quality control materials or serum pools that reasonably represent performance with patient specimens. Controls or serum pools should be used to monitor that the reagents are functioning properly and that correct procedures are being followed. An acceptable range for each lot of control material should be established by the laboratory. If control values are not within the expected range, follow normal trouble-shooting procedures. If assistance in the US is required, please call Sekisui Diagnostics Technical Marketing at (800) 332-1042. Outside the US, please contact your local distributor.

Quality control should be established in line with local, state and/or federal regulations and/or national regulations, or accreditation requirements.

## RESULTS

### Limitations / Interfering Substances

All interference studies were conducted according to CLSI EP7<sup>13</sup>.

Refer to the work of Young<sup>14</sup> for a review of the effects of drugs on clinical laboratory tests.

Heterophilic antibodies: Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed<sup>15</sup>.

### Expected Values

Samples from 196 apparently healthy individuals (101 males and 95 females) ranging in age from 18-79 were tested with the Sekisui Cystatin C assay using CLSI C28<sup>16</sup>. The results, calculated non-parametrically, were in the range 0.61-1.17 mg/L. This reference interval is consistent with literature values for cystatin C<sup>17</sup>.

Each laboratory should confirm the reference interval for the patient population it serves.


## SPECIFIC PERFORMANCE CHARACTERISTICS

### Accuracy

Comparative performance studies were conducted using the Sekisui Cystatin C on the Abbott Aeroset<sup>®</sup> clinical analyzer and the Dade Behring N Cystatin C method. 148 serum samples, with cystatin C concentrations between 0.48 and 7.73 mg/L were tested. The protocol followed the recommendations of CLSI EP9<sup>18</sup>.

The regression analysis is provided below:

Sekisui Cystatin C vs. Dade Behring N Cystatin C (n = 148)	
Slope	1.11
Intercept (mg/L)	0.099
Correlation Coefficient (r <sup>2</sup> )	0.997

### Precision

The precision of the Sekisui Cystatin C Reagent was determined using 3 levels of frozen, pooled human serum following CLSI EP5<sup>19</sup> by running samples in duplicate, twice per day, for 20 days on the Abbott Aeroset analyzer.

#### Within-Run Precision

Serum pool	Mean Recovery (mg/L)	Standard Deviation (mg/L)	CV
Level 1	0.72	0.01	1.7%
Level 2	2.79	0.02	0.7%
Level 3	5.06	0.08	1.5%

#### Total Precision

Serum pool	Mean Recovery (mg/L)	Standard Deviation (mg/L)	CV
Level 1	0.72	0.02	2.3%
Level 2	2.79	0.04	1.3%
Level 3	5.06	0.12	2.4%

### Assay Sensitivity

The Limit of the Blank (LOB) is the highest value 95% likely to be seen with a sample that contains no analyte. The LOB of the Sekisui Cystatin C Reagent on the Abbott Aeroset is 0.02 mg/L<sup>20</sup>.

The Limit of Detection (LOD) is the cystatin C concentration at which an observed test result is 95% likely to exceed the LOB. The LOD of the Sekisui Cystatin C Reagent on the Abbott Aeroset is 0.05 mg/L<sup>20</sup>.

The Limit of Quantitation (LOQ) represents the cystatin C concentration where the CV is  $\leq 10\%$ . The pooled standard deviation of the low end is 0.02mg/L, therefore the LOQ is 0.20mg/L.

### Specificity

The following substances, at the concentrations shown, did not affect the performance of the Sekisui Cystatin C assay in serum pools with cystatin C:

Theophylline	222 $\mu\text{mol/L}$
Tolbutamide	2.37 mmol/L
Trimethoprim	138 $\mu\text{mol/L}$

A high dose hook effect was not observed at concentrations of up to 25.0 mg/L cystatin C.

Rheumatoid factor was tested to 800 IU/mL and did not affect performance.

### Linearity


Using CLSI EP6<sup>21</sup>, the Sekisui Cystatin C method was linear from 0.2 to 8.0 mg/L. The samples were ad-mixtures of low and high serum. The high serum was prepared by supplementing a pool with purified cystatin C. Specimens with cystatin C concentrations above 8.0mg/L may be diluted at a ratio of 1:2 with physiological saline. Multiply the result by the dilution factor to obtain the cystatin C concentration in the sample.

<u>Substance</u>	<u>Concentration Tested</u>
Conjugated Bilirubin	50 mg/dL
Unconjugated Bilirubin	20 mg/dL
Hemoglobin	900 mg/dL
Triglyceride	940 mg/dL
Acetaminophen (paracetamol)	1324 $\mu\text{mol/L}$
Acetylsalicylic acid	3.26 mmol/L
Amoxicillin	206 $\mu\text{mol/L}$
Ampicillin	152 $\mu\text{mol/L}$
Atenolol	37.6 $\mu\text{mol/L}$
Azathioprine	10.8 $\mu\text{mol/L}$
Caffeine	308 $\mu\text{mol/L}$
Captopril	23 $\mu\text{mol/L}$
Chlorpheniramine maleate	0.8 mg/dL
Chlorpropamide	2.7 mmol/L
Cimetidine	79.2 $\mu\text{mol/L}$
Cyclosporin A	12 $\mu\text{mol/L}$
Digoxin	7.8 nmol/L
Doxycycline hyclate	67.5 $\mu\text{mol/L}$
Erythromycin	81.6 $\mu\text{mol/L}$
Furosemide	181 $\mu\text{mol/L}$
Glyburide	3.89 $\mu\text{mol/L}$
Ibuprofen	2425 $\mu\text{mol/L}$
Indomethacin	100 $\mu\text{mol/L}$
Loratidine	0.78 $\mu\text{mol/L}$
Losartan	16.6 mg/dL
Lovastatin	1.6 mg/dL
Metformin	310 $\mu\text{mol/L}$
Methotrexate	2 mmol/L
Methyl-DOPA	71 $\mu\text{mol/L}$
Metoprolol tartrate	18.7 $\mu\text{mol/L}$
Metronidazole	701 $\mu\text{mol/L}$
Nicotinic acid	2 mg/dL
Nitrofurantoin	16.8 $\mu\text{mol/L}$
Omeprazole	17.4 $\mu\text{mol/L}$
Phenytoin	198 $\mu\text{mol/L}$
Prednisone	0.84 $\mu\text{mol/L}$
Promethazine hydrochloride	4.22 $\mu\text{mol/L}$
Propranolol hydrochloride	7.71 $\mu\text{mol/L}$
Quinidine	37 $\mu\text{mol/L}$
Salicylic Acid	4.34 mmol/L
Simvastatin	0.8 mg/dL
Sirolimus	75 ng/mL
Sulfamethoxazole	1.58 mmol/L
Tacrolimus	50 nmol/L

#### The Americas

Sekisui Diagnostics, LLC  
31 New York Avenue  
Framingham, MA 01701 USA  
Phone: 800-332-1042  
Fax: 800-762-6311  
Email: [questions@sekisuidiagnostics.com](mailto:questions@sekisuidiagnostics.com)

#### International

Sekisui Diagnostics (UK), Limited  
50 Gibson Drive  
Kings Hill, West Malling  
KENT, ME19 4AF, UK  
Email: [info@sekisuidiagnostics.com](mailto:info@sekisuidiagnostics.com)

[www.sekisuidiagnostics.com](http://www.sekisuidiagnostics.com)

80-7138-00-02

5/12

©2012 Sekisui Diagnostics, LLC - All rights reserved.

Aeroset<sup>®</sup> is a registered trademark of Abbott Diagnostics.

**Definitions for Symbols / Erklärung der Symbole / Légende  
des symboles / Definizioni dei simboli / Significado de los  
símbolos / Definities van symbolen / Symboler /  
Symboldefinitioner**

This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.  
Dieses Produkt erfüllt die Anforderungen der Europäischen Richtlinie für In-vitro-Diagnostika.  
Ce produit remplit les conditions de la Directive européenne relative aux dispositifs médicaux de diagnostic *in vitro*.  
Questo prodotto soddisfa i requisiti della Direttiva Europea per i dispositivi medicali diagnostici in vitro  
Este producto cumple los requisitos de la Directiva Europea para dispositivos médicos de diagnóstico *in Vitro*  
Dit product voldoet aan de vereisten van de Europese Richtlijn voor in vitro diagnostische medische apparatuur.  
Denna produkt uppfyller kraven i det europeiska direktivet om medicintekniska produkter för in vitro-diagnostik.  
Produktet opfylder kravene i EU-direktivet for medicinsk udstyr til In Vitro diagnostik

Catalog number  
Katalognummer  
Réf. Catalogue  
Numero catalogo  
Número del catálogo  
Catalogusnummer  
Katalognummer  
Katalognummer

For *in vitro* diagnostic use  
Für die In-vitro-Diagnostik  
À usage diagnostique in vitro  
Per uso diagnostico in vitro  
Para uso diagnóstico *in Vitro*  
Voor in vitro diagnostisch gebruik  
För in vitro-diagnostik  
Til in vitro-diagnostisk brug

Temperature limitation  
Temperaturbegrenzung  
Limite de température  
Limiti di temperatura  
Limitación de temperatura  
Temperaturlægrenzing  
Temperaturgräns  
Temperaturgrænser

Manufactured by  
Hersteller  
Fabricant  
Prodotto da  
Fabricado por  
Geproduceerd door  
Tillverkad av  
Fremstillet af

Use by  
Verfallsdatum  
Date de péremption  
Utilizzare entro  
Fecha de caducidad  
Uiterste gebruiksdatum  
Sista användningsdag  
Anvendes inden

Batch code  
Chargennummer  
Code du lot  
Codice lotto  
Código del lote  
Batch-code  
Lotkod  
Batch-kode

Consult instructions for use  
Siehe Gebrauchsanweisung  
Consulter le mode d'emploi  
Consultare le istruzioni  
Consulte las instrucciones de uso  
Lees vóór gebruik de instructies documenten  
Se bruksanvisningen  
Læs brugsanvisningen

Caution, consult accompanying documents  
Achtung: Begleitdokumente lesen  
Attention : lire les documents joints  
Attenzione, consultare i documenti allegati per l'uso  
Atención, consulte la documentación adjunta  
Voorzichtig, lees de begeleverde  
Försiktighet – se bifogade dokument  
Forsigtig! se den medfølgende dokumentation

Authorized Representative in the European Community  
Autorisierter Vertreter in der Europäischen Gemeinschaft  
Représentant agréé en Europe  
Rappresentante autorizzato nella Comunità Europea  
Representante autorizado en la Comunidad Europea  
Geautoriseerd vertegenwoordiger in de Europese Unie  
Auktoriserad representant i EG  
Autoriseret repræsentant i EU

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1. Grubb, A. Lotberg H. Human gamma – trace, a basic microprotein; amino acid sequence and presence in the adenohypophysis. *Proc Natl Acad Sci USA* 1982; 79:3024-7.
2. Barrett AJ. Daviies ME, Grubb A. The place of human gamma-trace (Cystatin C) among the cysteine proteinase inhibitors. *Biochem Biophys Res Commun* 1984; 120: 631-6
3. Barrett AJ. Fritz H, Grubb A et al. Nomenclature and classification of the proteins homologous with chicken Cystatin. *Biochem J* 1986; 236:312.
4. Grubb A Diagnostic value of analysis of Cystatin C and protein HC in biological fluids. *Clin Nephrol* 1192;38: S20-7
5. Abrahamson M, Olafsson, I, Palsdttier A, et:al. Structure and expression of the human Cystatin C gene. *Biochem J* 1990; 268: 287-94.
6. Newman DJ. Thakkar H. Edwards RG, et al., Serum Cystatin C measured by automated immunoassay: a more sensitive marker of changes in GFR than serum creatinine. *Kidney Int* 1995; 47: 312-8.
7. Coll E, Botey A Alvarez L, et al Serum Cystatin C as a new marker for noninvasive estimation of glomerular filtration rate and as a marker for early renal impairment. *Am J Kidney Dis* 2000; 36: 29-34.
8. Leavering JH, Thai PJ, Van der Waart M, Schuurs AH. A sol particle agglutination assay for human chorionic gonadotrophin. *J Immunol Methods* 1981;45: 183-94.
9. Wielaard F. Denissen A. Van der Veen L, Rutjes I. A sol-particle immunoassay for determination of anti-rubella antibodies; development and clinical validation. *J. Virol Methods* 1987;17: 149-58.
10. Tanaka M, Matsuo K, Enomoto M, Mizuno K. A sol particle homogeneous immunoassay for measuring serum Cystatin C *Clin Biochem* 2004; 37: 27-35.
11. Richardson JH and Barkley WE, eds. *Biosafety in Microbiological and Biomedical laboratories*, U.S. Dept. of Health and Human Services, Public Health Service, HHS Publication No. (CDC) 84-8395, Washington, DC: 1984.
12. Clinical and Laboratories Standards Institute. *Procedures for the Handling and Processing of Blood Specimens: Approved Guideline*. CLSI document H18-A Villanova, PA: 1990.
13. Clinical and Laboratories Standards Institute. *Interference Testing in Clinical Chemistry: Approved Guideline*. CLSI document EP7-A2 Villanova, PA: 2005.
14. Young DS, *Effects of Drugs on Clinical Laboratory Tests*, 3<sup>rd</sup> ed. Washington, D.C.: AACC Press, 1990.
15. Boscato LM and Stuart MC, *Heterophilic Antibodies: A Problem for All Immunoassays*. *Clin Chem* 1988; 34: 27-33.
16. Clinical and Laboratories Standards Institute. *How to Define and Determine Reference Intervals in the Clinical Laboratory: Approved Guideline*. CLSI document C28-A2, Villanova, PA:2000.
17. Mussap, M and Plebani, M. *Critical Reviews in Clinical Laboratory Sciences*, 41(5-6):467-550, 2004.
18. Clinical and Laboratories Standards Institute. *Method Comparison and Bias Estimation Using Patient Samples: Approved Guideline*. CLSI document EP9-A2. Villanova, PA:2002.
19. Clinical and Laboratories Standards Institute. *Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline*. CLSI document EP5-A2. Villanova, PA:2004.
20. Clinical and Laboratories Standards Institute. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI document EP17-A. Villanova, PA:2004.
21. Clinical and Laboratories Standards Institute. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. NCCLS document EP6-A. Villanova, PA:2003.

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©2012 Sekisui Diagnostics, LLC - Alle rettigheder forbeholdt.