CHOLESTEROL-SL ASSAY

CATALOGUE NUMBER: 234-60 **SIZE:** 2 x 100 mL
234-99 1 x 1000 mL

INTENDED USE



For the IN VITRO quantitative measurement of total cholesterol in serum.

TEST SUMMARY

Cholesterol is a high molecular weight lipid found exclusively in animals. It is synthesized primarily in the liver and gut, and is released from those sites into the circulatory system⁽¹⁾. Serum cholesterol concentrations are affected by numerous factors, including hepatic and endocrine diseases, hereditary factors, and diet. Hypercholesterolemia is associated with increased risk of developing atherosclerosis.

Methodologies for measurement of cholesterol have been critically reviewed by Tonks⁽²⁾, as well as by Zak⁽³⁾. The enzymatic method used in this assay is a modification of that described by Allain et.al.⁽⁴⁾ and Roschlau et.al.⁽⁵⁾.

TEST PRINCIPLE

The cholesterol esters are hydrolyzed to free cholesterol by cholesterol esterase (CE). The free cholesterol is then oxidized by cholesterol oxidase (CO) to cholest-4-ene-3-one with the simultaneous production of hydrogen peroxide. The hydrogen peroxide produced couples with 4-aminoantipyrine and p-hydroxybenzoate, in the presence of peroxidase, to yield a chromogen with maximum absorbance at 505 nm.

The intensity of the color produced is directly proportional to the concentration of total cholesterol in the sample.

REAGENTS

Cholesterol Reagent: A buffered solution containing 0.5 mmol/L 4-aminoantipyrine, 25 mmol/L p-hydroxybenzoic acid, >240 U/L cholesterol esterase (mammalian and microbial), >150 U/L cholesterol oxidase (microbial), >1600 U/L peroxidase (botanical), surfactants, stabilizers, and preservatives.

WARNINGS & PRECAUTIONS FOR USE

Avoid contact with skin and eyes.
Harmful if swallowed.
See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use. Supplied reagent is stable at 2-8°C until expiry date. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum is the specimen of choice. The specimen should be drawn following a 12 hour fast.

SAMPLE STORAGE

Samples for the cholesterol assay are stable for one week at 2-8°C and for at least 4 weeks at -20°C.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁶⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia, and hemolysis were evaluated for this method on a Roche/Hitachi® 704 analyzer.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
161 mg/dL	4.1 mmol/L	Hemoglobin	500 mg/dL	78 µmol/L
184 mg/dL	4.8 mmol/L	Bilirubin	20 mg/dL	342 µmol/L
168 mg/dL	4.3 mmol/L	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

Interferences in the measurement of cholesterol using cholesterol oxidase procedures have been studied by Pesce and Boudourian⁽⁷⁾. A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁸⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' Cholesterol-SL reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Spectrophotometer or analyzer capable of accurately measuring absorbance at 505 nm.
2. Calibration material.
3. Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:100 and wavelength readings of 505 nm and 660 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the cholesterol concentration of each sample.

TEST LIMITATIONS

A sample with a cholesterol concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽⁹⁾

Values for selecting adults at moderate and high risk of coronary artery disease who require treatment.*

	Desirable	Borderline Risk for CHD	High Risk for CHD
Cholesterol Levels (Adults over age 20)	< 200 mg/dL (5.2 mmol/L)	200-239 mg/dL (5.2-6.2 mmol/L)	≥ 240 mg/dL (6.2 mmol/L)

*National Cholesterol Education Program (NCEP) Adult Treatment Panel Guidelines

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 704 analyzer unless otherwise stated.

RESULTS

Cholesterol concentration is reported as mg/dL (mmol/L).

REPORTABLE RANGE (CLSI EP6)⁽⁶⁾

The linearity of the procedure described is 600 mg/dL (15.5 mmol/L). The lower limit of detection of the procedure described is 1.2 mg/dL (0.03 mmol/L). This data results in a reportable range of 1.2 to 600 mg/dL (0.03 to 15.5 mmol/L).

PRECISION STUDIES (CLSI EP5)⁽⁶⁾

Data was collected on two concentrations of a control sera using a single lot of reagent in forty runs conducted over twenty days.

Concentration		Total SD		Total CV%	Within Run SD		Within Run CV%
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mg/dL	mmol/L
191.5	4.95	2.1	0.06	1.1	1.1	0.03	0.6
100.2	2.59	1.2	0.03	1.2	0.6	0.02	0.6

ACCURACY (CLSI EP9)⁽⁶⁾

The performance of this method (y) was compared with the performance of a similar cholesterol method (x) on a Roche/Hitachi® 704 analyzer. Forty-nine patient serum samples ranging from 132 to 403 mg/dL (3.4 to 10.4 mmol/L) were tested and gave a correlation coefficient of 0.9971. Linear regression analysis gave the following equation:

This method = 1.05 (reference method) – 13.9 mg/dL (0.36 mmol/L).

TRADEMARKS

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

The Americas
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

International
Sekisui Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK

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peidiagnosticttechnical@sekisuidiagnostics.com

Internacional

Sekisui Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, RU

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Temperature limitation
Límites de temperatura

REFERENCES/ REFERENCIAS

1. Burtis, C. and Ashwood, E., Editors, *Tietz Textbook of Clinical Chemistry*, W.B. Saunders Co., Toronto, 1004-1008 (1994).
2. Tonks, D.B., *The Estimation of Cholesterol in Serum, A Classification and Critical Review of Methods*, Clin. Biochem. 1, 12 (1967).
3. Zak, B., *Cholesterol Methodologies: A Review*, Clin. Chem. 23, 1201 (1977).
4. Allain, C.C., Poon, L., Chan, S.G., Richmond, W., Fu, P., *Enzymatic Determination of Total Serum Cholesterol*, Clin. Chem. 20, 470 (1974).
5. Roschlau, P., Bernt, E., Gruber, W., *Enzymatische Bestimmung des Gesamt-Cholesterins in Serum*, Z. Klin. chem. Klin. Biochem. 12, 226. (1974).
6. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
7. Pesce, M.A., Bodourian, S.H., *Interference with the Enzymatic Measurement of Cholesterol in Serum by Use of Five Reagent Kits*, Clin. Chem. 23, 757-760 (1977).
8. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Third Edition, Washington, 1990.
9. Rafai, N., Warnick, R.G., (Ed) *Laboratory Measurement of Lipids, Lipoproteins and Apolipoproteins*, AACC Press, Washington, DC, 25, (1994).

Authorized Representative/ Representante autorizado:
Sekisui Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ
United Kingdom/ Reino Unido
Tel: +44 (0) 1622 607800
Fax: +44 (0) 1622 607801

IN23460-14
October 28, 2015

The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.