EN

## DIRECT BILIRUBIN-SL ASSAY

**CATALOGUE NUMBER:** 247-30 **SIZE:** R1: 3 x 100 mL, R2: 1 x 75 mL  
 247-50A R1: 1 x 1000 mL  
 247-50B R2: 1 x 300 mL

### INTENDED USE

For the IN VITRO quantitative measurement of direct bilirubin concentration in serum.

### TEST SUMMARY

Bilirubin is a product of hemoglobin breakdown by the reticuloendothelial system and exists in two forms. Unconjugated (indirect) bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and excreted.

Elevated serum bilirubin is indicative of impairment of the liver, excessive hemolysis, or obstruction of the biliary tract.

Bilirubin determination is generally based on the reaction of bilirubin with diazotized sulfanilic acid, as described by Ehrlich.<sup>(1)</sup> In this method, direct (conjugated) bilirubin couples with a diazonium salt in the presence of sulfamic acid to form the colored compound, azo-bilirubin.

### TEST PRINCIPLE

Direct Bilirubin + 2,4-dichlorophenyldiazonium salt → Azobilirubin

Direct bilirubin reacts with a diazonium salt forming azobilirubin. The increase in absorbance at 546 nm due to azobilirubin formation is proportional to the direct bilirubin concentration.

### REAGENTS

DBili-SL Reagent (R1): A solution containing buffer (pH 1.20 at 25°C), 0.1 mmol/L EDTA-Na<sub>2</sub>, 154 mmol/L NaCl, and 100 mmol/L Sulfamic acid.

DBili-SL Diazo Reagent (R2): A solution containing buffer (pH 0.50 at 25°C), 0.5 mmol/L 2,4 dichlorophenyldiazonium salt, 0.9 mol/L HCl, and 0.13 mmol/L EDTA-Na<sub>2</sub>.

### WARNINGS & PRECAUTIONS FOR USE

Danger

Contains: Sulfamic acid; Hydrochloric acid

H314 – Causes severe skin burns and eye damage.

Prevention – P260 – Do not breathe vapour/spray.

P264 – Wash thoroughly after handling. P280 – Wear protective gloves/protective clothing/eye protection/ face protection.

Response – P301 + P330 + P331 – IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303 + P361 + P353 – IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

P363 – Wash contaminated clothing before reuse.

P304 + P340 – IF INHALED: Remove person to fresh air and keep comfortable for breathing.

P310 – Immediately call a POISON CENTRE or doctor/physician.

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Storage – P405 – Store locked up.

Disposal – P501 – Dispose of contents/container in accordance with local/regional/national/international regulations.

Avoid ingestion.

Avoid contact with skin and eyes.

See Material Safety Data Sheet for additional information.

### REAGENT, PREPARATION, STORAGE & STABILITY

Reagents are ready for use.

Supplied reagents are stable at 2-8°C until expiry date.

Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum. Samples should be protected from light as direct (conjugated) bilirubin is unstable.<sup>(1)</sup>

### SAMPLE STORAGE

Samples may be stored at 2-8°C for up to 3 days or at minus 70°C for up to 3 months in the dark.<sup>(1)</sup>

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(2)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interference from hemolysis and lipemia were evaluated for this direct bilirubin method on a Roche/Hitachi® analyzer using a significance criterion of >10% variance from control.

Concentration		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
mg/dL	µmol/L			
0.10	1.7	Intralipid	400 mg/dL	1200 mg/dL (13.5 mmol/L) Simulated Triglycerides

Hemoglobin produces significant interference with this method; therefore, hemolyzed samples should not be used.

Samples containing indocyanine green should not be used.

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(3)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Sekisui Diagnostics' Direct Bilirubin reagents.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material.
3. Quality Control materials.

#### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:16:4 and a wavelength reading of 546 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependant on the system and the parameters used.

#### QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state, and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

## CALCULATIONS

The analyzer automatically calculates the direct bilirubin concentration of each sample.

## TEST LIMITATIONS

A sample with direct bilirubin concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

## REFERENCE INTERVALS<sup>(1)</sup>

0.0-0.2 mg/dL (0.0-3.4 µmol/L) (adults and infants after 1 month).

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on an automated analyzer unless otherwise stated.

## RESULTS

Direct Bilirubin concentration is reported as mg/dL (µmol/L).

## REPORTABLE RANGE (CLSI EP6)<sup>(2)</sup>

The linearity of the procedure described is 10.0 mg/dL (171 µmol/L). The lower limit of detection of the procedure described is 0 mg/dL (0 µmol/L). This data results in a reportable range of 0-10 mg/dL (0-171 µmol/L).

## PRECISION STUDIES (CLSI EP5)<sup>(2)</sup>

Data was collected on two sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV %	Within Run SD		Within Run CV %
mg/dL	µmol/L	mg/dL	µmol/L		mg/dL	µmol/L	
0.25	4.3	0.018	0.31	7.3	0.005	0.09	2.1
3.22	55.1	0.139	2.38	4.3	0.014	0.24	0.4

## ACCURACY (CLSI EP9)<sup>(2)</sup>

The performance of this method (y) was compared with the performance of a similar method (x) on a Roche/Hitachi® analyzer. Fifty patient serum samples ranging from 0.06-10.8 mg/dL (1-185 µmol/L) gave a correlation coefficient of 0.9987. Linear regression analysis gave the following equation:

$$\text{This method} = 0.89 (\text{reference method}) + 0.132 \text{ mg/dL} (2.3 \text{ µmol/L}).$$

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**SEKISUI**  
DIAGNOSTICS

**The Americas**  
Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265  
Fax: 902-628-6504  
Email: questions@sekisuidiagnostics.com  
peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

**International**  
Sekisui Diagnostics (UK) Limited  
Liphook Way  
Allington, Maidstone  
KENT, ME16 0LQ, UK

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Temperature limitation  
Límites de temperatura

## REFERENCES/ BIBLIOGRAFÍA

1. Burtis, A. Carl, Ashwood, Edward R. (Editors) Tietz Textbook of Clinical Chemistry, Second Edition, W.B. Saunders Co., Toronto, (1994).
2. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
3. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Third Edition, Washington, 1990.

Authorized Representative/ Representante autorizado:  
Sekisui Diagnostics (UK) Limited  
Liphook Way  
Allington, Maidstone  
KENT, ME16 0LQ  
United Kingdom/ Reino Unido

The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.

IN24710-17  
November 1, 2016