LDL Direct Liquid Select Cholesterol Reagent

CATALOGUE NUMBER: 7120 **SIZE:** R1 1 x 30 mL
R2 1 x 10 mL

INTENDED USE

For the direct, quantitative measurement of low-density lipoprotein cholesterol (LDL-C) concentration in human serum or plasma.

TEST SUMMARY

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. The phospholipid, free cholesterol and protein constitute the outer surface of the lipoprotein particle, while the inner core contains mostly esterified cholesterol and triglyceride. These particles serve to solubilize and transport cholesterol and triglyceride in the bloodstream.

The relative proportions of protein and lipid determine the density of these lipoproteins and provide a basis on which to begin their classification.¹ These classes are: chylomicrons, very-low density lipoprotein (VLDL), low-density lipoprotein (LDL) and high density lipoprotein (HDL). Numerous clinical studies have shown that the different lipoprotein classes have very distinct and varied effects on coronary heart disease risk.²⁻⁴ The studies all point to LDL cholesterol as the key factor in the pathogenesis of atherosclerosis and coronary artery disease (CAD)²⁻⁸, while HDL cholesterol has been observed to have a protective effect. Even within the normal range of total cholesterol concentrations, an increase in LDL cholesterol can occur with an associated increased risk for CAD.⁴

TEST PRINCIPLE

The LDL Direct Liquid Select Cholesterol assay is a homogeneous method for directly measuring LDL-C concentrations in serum or plasma, without the need for any off-line pretreatment or centrifugation steps.

The method is in a two reagent format and depends on the properties of a unique detergent. This detergent (Reagent 1) solubilizes only the non LDL lipoprotein particles. The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non color forming reaction. A second detergent (Reagent 2) solubilizes the remaining LDL particles and a chromogenic coupler allows for color formation. The enzyme reaction with LDL-C in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

REAGENTS

Composition of Reagents:

Component	Ingredients	Concentration
Reagent 1	Buffer Detergent 1 Cholesterol esterase (Pseudomonas sp.) Cholesterol oxidase (Cellulomonas sp.) Peroxidase (Horseradish) 4-Aminoantipyrine Ascorbic oxidase (Curcubita sp.) Preservative	<1.0% <1500 U/L <1500 U/L <1300 ppg U/L <0.1% <3000 U/L
Reagent 2	Buffer (pH 6.3) Detergent 2 N,N-bis (4-sulfobutyl) -m-toluidine, disodium (DSBmT) Preservative	<1.0% <1.0 mM

WARNINGS AND PRECAUTIONS FOR USE

1. For *in vitro* diagnostic use.
2. Do not pipette by mouth.
3. All specimens used in this test should be considered potentially infectious. Universal precautions as they apply to your facility should be used for handling and disposal of materials during and after testing.
4. Do not use reagents after the expiration date printed on the label.

See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagent 1: Reagent 1 is ready to use as packaged.

Reagent 2: Reagent 2 is ready to use as packaged.

All unopened reagents are stable until the expiration date on the label when stored at 2-8°C.

Once opened, Reagent 1 and Reagent 2 is stable for 4 weeks at 2-8°C.

Once opened Reagent 1 and Reagent 2 have 4 weeks on board stability when stored at 2-8°C.

DO NOT FREEZE

REAGENT DETERIORATION

The following would indicate deterioration:

Inability to recover control values.
Presence of turbidity.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Patients are not required to fast prior to blood collection. Serum, EDTA-treated or heparinized plasma are the recommended specimens.

If not analyzed promptly, specimens may be stored at 2-8°C for up to 5 days. If specimens need to be stored for longer than 5 days, they may be stored frozen at -80°C.

ANALYTICAL SPECIFICITY

All interference studies were conducted according to NCCLS guideline No. EP7 for interference testing in clinical chemistry.¹²

Substance Tested Concentration with no significant (±10%) interference

Ascorbic Acid	50 mg/dL
Hemoglobin	500 mg/dL
Bilirubin	20 mg/dL
Gamma-Globulins	5000 mg/dL

1. Refer to the work of Young for a review of drug interference on serum LDL cholesterol levels.¹³
2. Protect the reagent from direct sunlight.
3. Anticoagulants containing citrate should not be used.
4. Patient samples should only be frozen once.
5. Samples with triglyceride values up to 1,293 mg/dL did not interfere with the results of the LDL Direct Liquid Select Cholesterol assay. Samples with Triglyceride levels >1,293 mg/dL should not be diluted.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Reagents and materials may be provided as follows:

Description	Configuration	Catalog Number
LDL Direct Liquid Select Cholesterol Reagent	R1 1 x 30 mL	7120
	R2 1 x 10 mL	


MATERIALS REQUIRED (BUT NOT PROVIDED)

Description	Configuration	Catalog Number
N-geneous® LDL Cholesterol Calibrator	3 x 1 mL	80-4610-02

1. Class A volumetric pipettes.
2. Distilled, deionized, Type II water or equivalent.
3. Analyzer capable of running two reagent chemistries.

TEST CONDITION

Below is a general example of the LDL Direct Liquid Select Cholesterol test procedure for a two reagent automated analyzer. All analyzer applications should be validated in accordance with NCEP and CLIA recommendations.^{10,11} For assistance with applications on automated analyzers, please contact Sekisui Diagnostics Technical Services at 800-565-0265.

CALIBRATION

The N-geneous® LDL Cholesterol Calibrator Kit is required for the calibration of this assay. Other commercially available LDL calibrators have not been tested with this assay and cannot be supported by Sekisui Diagnostics. The value of the LDL calibrator was assigned by procedures traceable to the CDC LDL cholesterol reference method.¹⁵ Refer to the instrument operator's manual for analyzer-specific calibration procedures and for guidance in determining calibration frequency.

Quality Control values should be within the expected range.

QUALITY CONTROL

Reliability of test results should be routinely monitored with quality-control materials or serum that reasonably represent performance with patient specimens.¹⁰ Controls or serum pools should be run with each assay to ensure that the reagents are functioning properly and that correct procedures have been followed. An acceptable range for each lot of control material should be established by the laboratory. If control values are not within the expected range, confirm procedures were performed correctly and follow normal troubleshooting measures. If assistance is required, call Sekisui Diagnostics Technical Services at 800-565-0265.

Quality control requirements should be established in accordance with local, state and/or federal regulations or accreditation requirements.

REFERENCE INTERVALS

The following NCEP cutpoints for patient classification are used for the prevention and management of coronary heart disease.⁸

It is recommended that each laboratory verify the reference interval for its patient population.

LDL Cholesterol	Classification
<130 mg/dL (<3.36 mmol/L)	Desirable
130-159 mg/dL (3.36 - 4.11 mmol/L)	Borderline High Risk
≥160 mg/dL (≥4.14 mmol/L)	High Risk

PERFORMANCE CHARACTERISTICS

RESULTS

To convert from conventional units to S.I. units, multiply the conventional units by 0.02586.¹⁴

$$\text{mg/dL} \times 0.02586 = \text{mmol/L LDL-cholesterol}$$

REPORTABLE RANGE

The limit of detection of the LDL Direct Liquid Select Cholesterol assay, quantified as 2 SDs plus the mean of twenty replicate measurements of saline, is 0.278 mg/dL on the Hitachi 911.

The LDL Direct Liquid Select Cholesterol test method is linear from 6.6 mg/dL to 992 mg/dL.

PRECISION STUDIES

Within-run precision for the LDL Direct Liquid Select Cholesterol Reagent was determined using three levels of frozen pooled human serum. Each run consisted of twenty replicate samples. Within-run precision studies produced the following results on the Hitachi 911:

Serum Pool	LOW	MID	HIGH
n	20	20	20
Mean LDL Cholesterol (mg/dL)	98.1	146.5	209.8
Standard Deviation (mg/dL)	0.72	0.96	1.31
Coefficient of Variation (%)	0.73	0.66	0.62

Between-run precision was determined using three levels of frozen pooled human serum. The LDL Direct Liquid Select Cholesterol assay was run twice per day in duplicate over 10 days. Between-run precision studies produced the following results on the Hitachi 911:

Serum Pool	LOW	MID	HIGH
n	40	40	40
Mean LDL Cholesterol (mg/dL)	98.1	142.7	207.3
Standard Deviation (mg/dL)	2.2	2.8	3.6
Coefficient of Variation (%)	2.27	1.95	1.73

ACCURACY

Accuracy of LDL Direct Liquid Select Cholesterol Reagent method was verified by comparison to the Reference Method (Ultracentrifugation and cholesterol analysis)¹⁰ and the Direct LDL Immunoseparation method.

Studies comparing the LDL Direct Liquid Select Cholesterol method to the Reference Method (Ultracentrifugation) produced the following results:

Method	LDL Direct Liquid Select Cholesterol	Reference Method
n	54	54
Mean LDL Cholesterol (mg/dL)	122.5	125.1
Standard Deviation (mg/dL)	30.7	30.9
Regression Analysis	$y = 0.95x + 3.02 \text{ mg/dL}$	
Correlation Coefficient	$r = 0.96$	

Studies comparing the LDL Direct Liquid Select Cholesterol method to the Direct LDL cholesterol Immunoseparation method produced the following results:

Method	LDL Direct Liquid Select Cholesterol	Direct LDL
n	92	92
Mean LDL Cholesterol (mg/dL)	120.0	122.8
Standard Deviation (mg/dL)	30.5	31.6
Regression Analysis	$y = 0.94x + 4.46 \text{ mg/dL}$	
Correlation Coefficient	$r = 0.97$	

All trademarks, brands, product names and trade names are the property of their respective companies.

The Americas

Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265

Fax: 902-628-6504

Email: questions@sekisuidiagnostics.com

peidiagnostictechnical@sekisuidiagnostics.com

International

Sekisui Diagnostics (UK) Limited

Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK

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In vitro diagnostic medical device
Appareil médical de diagnostic *in vitro*.

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Caution, consult accompanying document
Prudence, consulter les documents joints

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1. Gotto AM, Lipoprotein metabolism and the etiology of hyperlipidemia, Hospital Practice, 23: Suppl.1, 4 (1988).
2. Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26; 566 (1985).
3. Badimon JJ, Badimon L., Fuester V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85:1234-41 (1990).

4. Castelli WP et al., Cholesterol and other lipids in coronary heart disease, Circulation, 55; 767 (1977).
5. Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med.,11;480 (1951).
6. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am.J. Med.,62;707 (1977).
7. William P., Robinson D., Baily A., High density lipoprotein and coronary risk factor, Lancet, 1;72 (1979).
8. Kannel WB, Castelli WP, Gordon T., Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Am. Intern. Med., 90; 85 (1979).
9. National Institutes on Health publication No. 93-3095, September 1993.
10. Bachorik PS et al, National Cholesterol Education Program Recommendations for Measurement of Low- Density Lipoprotein Cholesterol; Executive Summary, Clinical Chemistry, Vol. 41, No. 10, (1995).
11. Grundy SM et al., Summary of the Second Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II), JAMA, 269; 23; 3015-3023 (1993).
12. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation Protocol Number 7, Vol. 4, No.8, June 1984.
13. Young DS, Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACC Press, Washington DC, 3-104 thru 3-106, (1990).
14. Tietz NW, Clinical Guide to Laboratory Tests, W.B. Saunders Co., Philadelphia, p. 256, (1986).
15. National Reference System for Cholesterol. CRMLN LDL Cholesterol Protocol, May 2004.

N-geneous® is a registered trademark of Sekisui Diagnostics, LLC.
©2011 Sekisui Diagnostics, LLC - All rights reserved.

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The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.

IN7120-1
April 20, 2015

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N-geneous® est une marque de commerce déposée de Sekisui Diagnostics, LLC
©2011 Sekisui Diagnostics, LLC - Tous droits réservés.