

# osom<sup>®</sup>

## ImmunoDip<sup>®</sup>

### Urinary Albumin Test

TO AID IN THE EARLY DETECTION OF KIDNEY DISEASE




FOR USE WITH URINE

FEATURES:



- ✓ **ACCURATE;** Sensitivity > 99%; Specificity 98%
- ✓ **FAST;** Results in 3 minutes
- ✓ **EARLY DETECTION;** Sensitive at low levels of microalbumin (18 mg/L)
- ✓ **EASY-TO-USE;** Simply place test device into urine sample
- ✓ **ROOM TEMPERATURE STORAGE;** 18-month shelf life

TEST PROCEDURE / RESULTS



ORDERING INFORMATION

CONFIGURATION

CATALOG NUMBER

○ OSOM<sup>®</sup> IMMUNODIP<sup>®</sup> URINARY ALBUMIN TEST  
 ○ OSOM<sup>®</sup> IMMUNODIP<sup>®</sup> CONTROL KIT

25 TESTS

790-25 (CE MARKED)  
 791 (CE MARKED)



**THE AMERICAS**  
 Sekisui Diagnostics, LLC  
 31 New York Avenue  
 Framingham, MA 01701  
 Phone: 800 332 1042  
 Fax: 800 762 6311  
 Email: questions@sekisuidiagnostics.com  
 ISO 13485:2003 Certified

**INTERNATIONAL**  
 Sekisui Diagnostics (UK) Limited  
 50 Gibson Drive  
 Kings Hill, West Malling, Kent, ME19 4AF, UK  
 Phone: +44 (0) 1732 220022  
 Fax: +44 (0) 1732 220024  
 Email: info@sekisuidiagnostics.com  
 ISO 9001:2000 Certified, ISO 13485:2003 Certified



OSOM<sup>®</sup>, ImmunoDip<sup>®</sup> and the OSOM Training logo are registered U.S. trademarks of Sekisui Diagnostics, LLC.  
 ©2011 Sekisui Diagnostics, LLC – All rights reserved.