

# OSOM<sup>®</sup>

## Influenza A&B Test



FOR USE WITH NASAL SWAB SPECIMENS


FEATURES:

- ✓ **ACCURATE, DIFFERENTIAL DIAGNOSIS;** Differentiates between influenza A&B
- ✓ **QC INSIDE;** Includes 2 additional test sticks for external QC testing
- ✓ **NASAL SWAB SPECIMENS**
- ✓ **FAST;** Less than 1 minute hands on time; results in 10 minutes or less
- ✓ **ROOM TEMPERATURE STORAGE**
- ✓ **DETECTS H1N1 STRAINS\***

\*Although this test has been shown to detect the Mexico/4108/2009 H1N1 virus cultured from a positive human respiratory specimen, the performance characteristics of this device with clinical specimens that are positive for the 2009 H1N1 influenza virus have not been established.



TEST PROCEDURE / RESULTS



ORDERING INFORMATION

CONFIGURATION

CATALOG NUMBER

○ OSOM<sup>®</sup> INFLUENZA A&B TEST  
INFLUENZA A&B CONTROL KIT

25 Tests

190 or 190E (CE marked)  
191 or 191E (CE marked)



**THE AMERICAS**  
31 New York Avenue  
Framingham, MA 01701-9322  
Phone: 800 332 1042  
Fax: 800 762 6311  
Email: usdiagcustomerservice@genzyme.com  
ISO 13485:2003 Certified

**INTERNATIONAL**  
50 Gibson Drive  
Kings Hill, West Malling  
Kent, ME19 4AF, UK  
Phone: +44 (0) 1732 220022  
Fax: +44 (0) 1732 220024/5  
Email: ukdiagcustomerservice@genzyme.com  
ISO 9001:2000 Certified, ISO 13485:2003 Certified

