



## SPECIFICATION

ISO: 75F41G06-85  
March 1, 2009  
[ 5<sup>th</sup> edition ]

Product Name: L- $\alpha$ -GLYCEROPHOSPHATE OXIDASE [GPOM]  
(L-Glycero-3-phosphate:oxygen oxidoreductase, EC 1. 1. 3. 21)

Product Number: T-107

Origin: *Streptococcus* sp.

Description: Yellowish lyophilized powder

Activity:  $\geq 15.0$  U/mg

Appearance of solution: Clear  
(50U/mL in buffer solution)  
A400nm  $\leq 0.150$   
A500nm  $\leq 0.060$   
A650nm  $\leq 0.015$

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.*  
*Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE  
Plant Director  
Ohito diagnostics plant

H. Nagai

H.NAGAI  
Quality control Assistant Manager  
Ohito diagnostics plant