



Diagnostics Department
1-105 Kanda Jinbocho, Chiyoda-ku, Tokyo 101-8101, Japan
Phone +81-(0)3-3296-3618, Fax +81-(0)3-3296-3682
diagnostics @om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41L09-85
March 1, 2009
[4th edition]

Product Name: **LIPASE [LPM]**
(Triacyl glycerol acylhydrolase, EC 3. 1. 1. 3)

Product Number: **T-116**

Origin: *Chromobacterium viscosum*

Description: White to pale brownish lyophilized powder

Activity: 1500~3500 U/mg

Appearance of solution: Clear
(1mg/mL in buffer solution)
A400nm ≤0.100
A500nm ≤0.050
A650nm ≤0.050

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant