



*FAST / SIMPLE / ACCURATE*

# The FastPack® IP Free PSA Immunoassay

The FastPack® IP Free PSA Immunoassay, for use on the FastPack® IP System, allows any testing facility to produce point-of-care results in 12 minutes with less than one minute of hands-on time.

- Enhanced quality of service for patients
- Ability to monitor and treat in one visit
- Improved patient compliance and outcome
- Laboratory accuracy and precision



**SEKISUI**  
DIAGNOSTICS

*Because every result matters™*

# Simple Process



**1**  
*Dispense sample into the FastPack® Injection Port*



**2**  
*Place FastPack® into the FastPack® IP Instrument*



**3**  
*Press the START button*

## The FastPack® IP Free PSA Assay Specifications:

- Quantitative measurement of Free PSA
- Sample: serum
- Lower limit: 0.05 ng/mL
- Upper limit: 20 ng/mL
- Assay time: 12 minutes

ORDERING INFORMATION	CONFIGURATION	CATALOG NUMBER
FastPack® IP Free Prostate Specific Antigen (Free PSA)	30 Tests	25000045
FastPack® IP Free PSA Calibration	Kit	25000008
FastPack® PSA/Free PSA Control	Kit	25000003

### THE AMERICAS

Sekisui Diagnostics, LLC  
One Wall Street  
Burlington, MA 01803  
Phone: 800 332 1042  
Fax: 800 762 6311  
Email: [questions@sekisui-dx.com](mailto:questions@sekisui-dx.com)

### INTERNATIONAL

Sekisui Diagnostics (UK) Limited  
Liphook Way, Allington Maidstone,  
Kent ME16 0LQ, UK  
Phone: +44 1622 607800  
Fax: +44 1622 607801  
Email: [info@sekisui-dx.com](mailto:info@sekisui-dx.com)

[www.sekisuidiagnostics.com](http://www.sekisuidiagnostics.com)

**SEKISUI**  
DIAGNOSTICS

*Because every result matters™*