

SILARIS®

RSV

This single test **accurately aids in the diagnosis of RSV** using the trusted performance of **PCR technology**.



RSV Facts

ALMOST ALL CHILDREN WILL BE INFECTED WITH RSV BY THEIR SECOND BIRTHDAY.

RSV USUALLY CAUSES MILD RESPIRATORY ILLNESS, BUT IT CAN BE SERIOUS, ESPECIALLY IN INFANTS AND OLDER ADULTS.


IT IS THE MOST COMMON CAUSE OF BRONCHIOLITIS AND PNEUMONIA IN CHILDREN YOUNGER THAN 1 YEAR OF AGE.¹

Test





Simple 3-step process fits easily into your existing workflow



accurate

- ✓ Molecular results for rapid and effective patient management
- ✓ Accurate diagnosis helps facilitate proper infection control measures
- ✓ 90.2% sensitivity/95.6% specificity

affordable

- ✓ Single test for diagnosis
- ✓ Controls in every kit
- ✓ Molecular reimbursement (CPT Code 87634QW)

simple

- ✓ Nasal swab collection
- ✓ Compact and portable to facilitate near patient testing
- ✓ No age limitations
- ✓ 30 minutes to test results

ORDERING INFORMATION	Part No.
Silaris RSV Test	1030
Silaris Dock	1026
Silaris RSV Control Kit	1029

THE AMERICAS
 Sekisui Diagnostics, LLC
 One Wall Street
 Burlington, MA 01803
 Phone: 800 332 1042
 Fax: 800 762 6311
 Email: questions@sekisui-dx.com

INTERNATIONAL
 Sekisui Diagnostics (UK) Limited
 Liphook Way, Allington
 Maidstone, Kent, ME16 0LQ, UK
 Phone: +44 1622 607800
 Fax: +44 1622 607801
 Email: info@sekisui-dx.com



Because every result matters™

www.sekisuidiagnostics.com