

Stronger Clinical Performance Takes Lateral Flow Testing To The Next Level

osom ULTRA PLUS FLU A&B Test



Providing superior rapid results at the point-of-care. Fast, easy, cost effective so you can test and treat in one visit.



High Performance

Equivalent to or exceeding the performance of reader devices, without the need for an instrument



Results in 10 minutes



OSOM® Custom Care

Exceptional Support/Training by licensed medical technologists and experienced healthcare professionals




Made in the USA


Award-winning supply chain capabilities



SENSITIVITY



SPECIFICITY



SIMPLE TO USE / RESULTS IN MINUTES

osom[®]

ULTRA PLUS
FLU A&B Test

1 INSERT SWAB IN
BUFFER AND ROTATE

10X

2 INSERT TEST STICK
INTO BUFFER

3 READ RESULTS
AT 10 MINUTES



Two extra test sticks provided at no additional cost for external QC!

INTERPRETATION OF RESULTS

INFLUENZA A+



INFLUENZA B+



NEGATIVE



INVALID



#1032 - OSOM Ultra Plus Flu A&B Test Kit – 25 tests per kit

#1034 - OSOM Ultra Plus Flu A&B Control Kit – 5 Flu A Swabs, 5 Flu B Swabs

THE AMERICAS

Sekisui Diagnostics, LLC
One Wall Street
Burlington, MA 01803
Phone: 800 332 1042
Fax: 800 762 6311
Email: questions@sekisui-dx.com

INTERNATIONAL

Sekisui Diagnostics (UK) Limited
Liphook Way, Allington
Maidstone, Kent, ME16 0LQ, UK
Phone: +44 1622 607800
Fax: +44 1622 607801
Email: info@sekisui-dx.com

SEKISUI
DIAGNOSTICS

Because every result matters™

www.sekisuidiagnostics.com