



Total Protein measurements are used as an aid to monitor and diagnose **liver** and **renal** diseases.

METHOD: BIURET; ENDPOINT

FEATURES:

FOR THE QUANTITATIVE MEASUREMENT OF TOTAL PROTEIN

Total Protein

The Sekisui Total Protein assay uses a Biuret procedure as a reliable determinant for the quantitative measurement of total protein in serum.

- ✓ No significant interference displayed from samples with elevated levels of icterus
- ✓ One part, stable liquid
- ✓ Applicable to multiple instrument platforms

BENEFITS:

- ✓ Reduces the need for sample dilutions
- ✓ Easy to use, no additional preparation required
- ✓ Laboratory flexibility

ORDERING INFORMATION

CONFIGURATION

CATALOG NUMBER

TOTAL PROTEIN

R1 4 X 125 mL

200-55

DC-CAL CALIBRATOR

5 X 3 mL

SE-035

DC-TROL LEVEL 1

10 X 5 mL

SM-052

DC-TROL LEVEL 2

10 X 5 mL

SM-056



THE AMERICAS
 Sekisui Diagnostics, LLC
 4 Hartwell Place
 Lexington, MA 02421
 Phone: 800 332 1042
 Fax: 800 762 6311
 Email: questions@sekisuidiagnostics.com

INTERNATIONAL
 Sekisui Diagnostics (UK) Limited
 50 Gibson Drive
 Kings Hill, West Malling, Kent, ME19 4AF, UK
 Phone: +44 (0) 1732 220022
 Fax: +44 (0) 1732 220024
 Email: info@sekisuidiagnostics.com

