



SPECIFICATIONS

ANALYSIS	
TEST THROUGHPUT	Up to 400 tests/hour (clotting assays); Up to 200 tests/hour (chromogenic and immunoturbidimetric assays)
SIMULTANEOUS MEASUREMENT ITEMS	Up to 20 items
MEASUREMENT PRINCIPLE	Light scattering photometry method for clotting assays; Absorbance photometry method for chromogenic and latex immunoturbidimetric assays
LIGHT SOURCE	LED (660nm), halogen lamp (405/570/730nm)
CALIBRATION	Automatic dilution setting; Multi-point calibration curve; INR calibration curve; Automatic switching of calibration curves
QUALITY CONTROL	Daily time-series, daily twin plot
SAMPLES	
SAMPLE LOADING	Continuous rack loading of up to 10 racks of 5 samples
STAT SAMPLES	2 dedicated STAT lanes
CAPABILITY TO LINK TO LAS	Point in space sampling — Option
CLOSED TUBE SAMPLING (CTS)	Cap piercing — Option
SAMPLE DISPENSING	Sample volumes from 2-100µl (0.5µl increments)
SAMPLE BARCODE READER	Internal reader, reads 16 digit, 2 of 5, CODE 39, CODE 128, NW 7 barcode types
REAGENTS	
REAGENT COMPARTMENT	30 refrigerated positions; 9 additional positions for deficient plasmas, diluent, and detergent
REAGENT BARCODE READER	Internal datamatrix code reader and handheld 2D code reader
OTHER	
OPERATING SYSTEM	Windows® 7
PHYSICAL DIMENSIONS /WEIGHT	Without CTS W 590 x D 740 x H 1190 mm / Approx. 185kg (Analysis unit only); With CTS W 590 x D 915 x H 1190mm / Approx. 200kg (Analysis unit only)
DATA	Data storage for up to 2.5 years



Because every result matters™

THE AMERICAS
Sekisui Diagnostics, LLC
4 Hartwell Place
Lexington, MA 02421
Phone: 800 332 1042
Fax: 800 762 6311
Email: questions@sekisui-dx.com
ISO 13485 Certified



INTERNATIONAL
Sekisui Diagnostics (UK) Limited
Liphook Way, Allington
Maidstone, Kent, ME16 0LQ
Phone: +44 (0) 1732 220022
Fax: +44 (0) 1732 220024
Email: info@sekisui-dx.com
ISO 9001 Certified, ISO 13485 Certified

NOT FOR SALE IN THE USA

Because every result matters™ is a trademark of Sekisui Diagnostics.
CP3000 is a mark of Sekisui Medical Co., Ltd.
© 2017 Sekisui Diagnostics, LLC- All rights reserved.