



SPECIFICATION

ISO: 75F41C14-85

March 1, 2009

[4th edition]

Product Name:	CHOLESTEROL OXIDASE [CON II -FD] (Cholesterol:oxygen oxidoreductase, EC 1. 1. 3. 6)
Product Number:	T-101
Origin:	<i>Microorganism</i>
Description:	Yellowish lyophilized powder
Activity:	≥ 15 U/mg
Appearance of solution: (5U/mL in buffer solution)	Clear A400nm ≤ 0.100 A500nm ≤ 0.050 A650nm ≤ 0.010
Contaminating enzymes:	
Catalase	≤ 0.5 % (U/U)
Cholesterol esterase	≤ 0.02 % (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M.TOBE
Plant Director
Ohito diagnostics plant

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant