



Diagnosics Department
1-105 Kanda Jinbocho, Chiyoda-ku, Tokyo 101-8101, Japan
Phone +81-(0)3-3296-3618, Fax +81-(0)3-3296-3682
diagnosics @om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41P04-85
March 1, 2009
[4th edition]

Product Name: **PHOSPHOLIPASE-C [PLC]**
(Phosphatidylcholine cholinephosphohydrolase, EC 3. 1. 4. 3)

Product Number: **T-11**

Origin: *Bacillus cereus*

Description: White to brownish lyophilized powder

Activity: ≥ 30 U/mg

Appearance of solution: Clear
(10U/mL in buffer solution)
A400nm ≤ 0.100
A500nm ≤ 0.050
A650nm ≤ 0.025

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant