



Diagnosics Department  
1-105 Kanda Jinbocho, Chiyoda-ku, Tokyo 101-8101, Japan  
Phone +81-(0)3-3296-3618, Fax +81-(0)3-3296-3682  
diagnostics @om.asahi-kasei.co.jp

## SPECIFICATION

ISO: 75F41M01-85  
March 1, 2009  
[ 4<sup>th</sup> edition ]

Product Name: **MONOGLYCERIDE LIPASE [MGLP II]**  
(Glycerol-monoester acylhydrolase, EC 3. 1. 1. 23)

Product Number: **T-117**

Origin: *Microorganism*

Description: White lyophilized powder

Activity:  $\geq 20$  U/mg

Appearance of solution: Clear  
(1mg/mL in buffer solution)  
A400nm  $\leq 0.015$   
A500nm  $\leq 0.015$   
A650nm  $\leq 0.015$

Contaminating enzymes:  
Catalase  $\leq 0.5$  % (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.*  
*Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

-----  
M.TOBE  
Plant Director  
Ohito diagnostics plant

H. Nagai

-----  
H.NAGAI  
Quality control Assistant Manager  
Ohito diagnostics plant