



Diagnosics Department  
1-105 Kanda Jinbocho, Chiyoda-ku, Tokyo 101-8101, Japan  
Phone +81-(0)3-3296-3618, Fax +81-(0)3-3296-3682  
diagnosics @om.asahi-kasei.co.jp

## SPECIFICATION

ISO: 75F41A02-85  
March 1, 2009  
[ 4<sup>th</sup> edition ]

Product Name: **ACYL-COA OXIDASE [ACOD]**  
(Acyl-CoA: oxygen oxidoreductase, EC 1. 3. 3. 6)

Product Number: **T-17**

Origin: *Arthrobacter* sp.

Description: Yellowish lyophilized powder

Activity: 20~40 U/mg

Appearance of solution: Clear  
(20U/mL in buffer solution)  
A400nm ≤0.080  
A500nm ≤0.040  
A650nm ≤0.020

Contaminating enzymes:  
Catalase ≤1.00 (U/U)  
Glucose oxidase ≤0.01 % (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.  
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE  
Plant Director  
Ohito diagnostics plant

H. Nagai

H.NAGAI  
Quality control Assistant Manager  
Ohito diagnostics plant